

Surgical Probe Cover with Pull-Up

MADE IN
USA




- Designed for ultrasound transducers used in the operating room
- Features pull-up system
- Telescopically folded
- Sterile, single-use
- Made of latex-free material
- 24/box
- Total length of 96"
- Made in the USA

Our surgical probe covers with pull-up are an ideal solution for protecting patients and medical staff in the operating room. The cover's pull-up system reduces the risk of cross-contamination by allowing easy insertion of the probe. Its length of 96" allows medical staff to cover both the transducer and the cable, thus providing a maximum level of protection by preserving the sterile field. The soft portion of the cover protects the transducer and is made of a soft and resistant polyurethane (latex-free) material. The clear portion of the cover protects the cable and is clear polyethylene (latex-free). By using this cover, the cable of the transducer is free to run over patients and other surfaces in the sterile field without compromising it.

Item number	Description	Pckg	Price
103.0008	Surgical probe covers with pull-up - 5" x 96"	24/box	\$225

Telescopic Surgical Probe Cover

- € Designed for ultrasound transducers used in the operating room
- € Kit includes covers, sterile gel, and elastic bands
- € Individually wrapped and sterilized with EO
- € Telescopically folded
- € Sterile, single-use
- € Made of latex-free polyethylene material
- € Total length of 96"
- € 24/box
- € Made in the USA

The telescopic folding of this surgical cover is ideal for the operating room as it allows medical staff to easily invert the sheath over their hand, grab the transducer, and simply unfurl the rest of the cover. This sterile cover also comes with elastic bands, perfect to ensure a snug fit of the cover over the probe. In addition, its length provides coverage for both the transducer and the cord, thereby preserving the sterile field.

Item number	Description	Pckg	Price
103.0007	Surgical probe covers with pull-up - 5" x 96"	24/box	\$215