

Latex probe covers



- General purpose ultrasound probe cover
- Sterile latex material provides a cost-effective solution
- Designed for ultrasound-guided interventional procedures
- Cover can accommodate both linear and curved probes
- Each kit includes two elastic bands and a Parker sterile gel packet
- Durable and tear-resistant design
- Cover measures 1.5" x 7.9" (3.5 x 20 cm)
- Made in the USA

Our latex probe covers provide a cost-effective solution for your patient safety needs. This tear-resistant and durable cover creates a sterile barrier between the probe and the equipment, thereby reducing the risk of cross-contamination. In addition, it can accommodate both linear and curved probes as it is stretchable up to a diameter of 3.5" (8.89 cm).

Item number	Description	Pckg	Price
101.0027	1.5" x 7.9" (3.5 x 20 cm) latex sterile covers	40/box	\$147.00