

QuadLife Accessories

Home Environmental Control - Q-HEC

The Q-HEC (Home Environmental Control) provides the QuadJoy mouse with: 2 External Switch inputs, 1 extra USB port (useful for XBOX joystick function), the IR blaster port to remote control most TV's, provide a Nurse call function, 2 relay contacts to control devices, (future) a R/C tether to team with Radio Controlled hobby devices and the ability to battery power your QuadJoy. The Joystick option is Free to US Veterans, Hospital's, and Rehabilitation facilities.



Mini DIN 8-pin Cable

- **6-Feet - QC-DIN8-6**
- **10-Feet - QC-DIN8-10**
- **25-Feet - QC-DIN8-25**

Mini DIN 8-pin Cable is needed for connecting the QuadJoy mouse to either the External Jack Interface or the Wheelchair Jack Interface.



XBOX Compatible Joystick - QS-XBOX

Please enter the unit number at checkout for speedy delivery. The location of the unit number is at the bottom of the mouse.



External Jack Interface - Q-JACK

Allows 3 External Switches that can trigger any Action. The External Jack interface connects to the QuadLife mouse with a QC-DIN8-xx cable ordered separately available in 6', 10' and 25'.



Wall Power Supply - Q-Power-10

Wall Power Supply is an inexpensive, portable, alternative power source for your QuadJoy. If you don't like constantly draining your computer's battery and lugging it around to power the QuadJoy, the Wall Power Supply may be a great option, especially if you use Bluetooth.

