

Blom-Singer®

StomaSoft™ Laryngectomy Tube

Introducing the new StomaSoft™
HME-compatible Laryngectomy Tube



Features include:

- Size 8, size 9, size 10 and size 12
- Available in two lengths (36mm and 55mm)
- Designed with a reinforced strap holder
- Compatible with InHealth Technologies' HME systems and accessories
- Available with or without fenestrations

Feel the flexibility and see the clarity with the
Blom Singer® **StomaSoft™** Laryngectomy Tube

Blom-Singer®

StomaSoft™ Laryngectomy Tube

SIZE	ID/OD	FENESTRATION	PART NUMBER
8/36	9.5/12.0	Non-fenestrated	BE 6698
8/55	9.5/12.0	Non-fenestrated	BE 6699
9/36	10.5/13.5	Non-fenestrated	BE 6700
9/55	10.5/13.5	Non-fenestrated	BE 6701
10/36	12.0/15.0	Non-fenestrated	BE 6702
10/55	12.0/15.0	Non-fenestrated	BE 6703
12/36	13.5/17.0	Non-fenestrated	BE 6704
12/55	13.5/17.0	Non-fenestrated	BE 6705
8/36	9.5/12.0	Fenestrated	BE 6698F
8/55	9.5/12.0	Fenestrated	BE 6699F
9/36	10.5/13.5	Fenestrated	BE 6700F
9/55	10.5/13.5	Fenestrated	BE 6701F
10/36	12.0/15.0	Fenestrated	BE 6702F
10/55	12.0/15.0	Fenestrated	BE 6703F
12/36	13.5/17.0	Fenestrated	BE 6704F
12/55	13.5/17.0	Fenestrated	BE 6705F

INHEALTH
TECHNOLOGIES®

For ordering information, go to
www.inhealth.com

©2019 InHealth Technologies – Manufactured by Freudenberg Medical, LLC (190512.00)