

XL Performance

Thanks to its modern and smart power management, the SERVOX® digital XL speech aid offers significantly extended battery performance.

XL Flexibility

SERVOX® digital XL provides more freedom and flexibility due to already included rechargeable batteries, a universal travel adapter and a high-quality belt pouch.

XL Quality

The SERVOX® digital XL speech aid combines technologies proven since 1958 with the state-of-the-art components and first-class raw materials.

Technical data: Battery Charger

| | |
|------------------|-----------------------------|
| VARTA Type | 57677 |
| Cell sizes | AAA, AA, 9V, Device |
| System to charge | Nickel Metal Hydride (NiMH) |
| Charging time | 4 hours** |

Technical data: Travel Adapter

| | |
|---------------|-------------------------------|
| Voltage | 100 V – 240 V |
| Maximum power | 660 W @ 100 V, 1350 W @ 230 V |

SERVOX® digital XL set includes

SERVOX® digital XL speech aid, Oral tube with adapter, Screw cap „soft tone“, Carrying strap with replacement strap protection, Belt pouch, AAA rechargeable batteries, Charger (commercially available), travel adapter, 2 battery holders, Instruction manual

Order Information

| Description | Unit | Order No. |
|---|------|-----------|
| Speech Aid SERVOX® digital XL | 1 | 37270 |
| Standard screw cap, completely assembled | 1 | 22774 |
| Screw cap "soft tone", completely assembled | 1 | 23285 |
| Oral tube with adapter | 1 | 18253 |
| Strap set (carrying strap) | 1 | 17093 |
| Belt pouch | 1 | 16555 |
| Bottom cover | 1 | 37505 |
| O-ring | 1 | 14831 |
| Battery holder | 1 | 37574 |
| Rechargeable batteries (AAA) | 1 | 37561 |
| Charger | 1 | 37469 |

* Average battery life under factory settings in continuous operation | ** Depends on cell capacity: Based on 4AA 1600 mAh



SERVOX® digital XL
Are you ready to speak?


SERVOX® digital XL

Are you ready to speak?

XL Performance – XL Flexibility – XL Quality

The **SERVOX® digital XL** speech aid is a high-quality new edition of the well-known **SERVOX® digital**. With its modern power management the battery performance has been significantly increased, extending speaking time by more than 50%* compared to the previous model. **SERVOX® digital XL** provides up-to-date vocal rehabilitation for laryngectomized and tracheotomized patients.

Battery life



* Average battery life under factory settings in continuous operation

Optimization of the battery performance




The new power management of the **SERVOX® digital XL** includes rechargeable batteries, which leads to the 50%*- increase in battery life compared to the previous model. The standard AAA-batteries create new form of freedom and flexibility.

Smart charging management



The new high-end charger of **SERVOX® digital XL** is easy to handle and has an integrated USB-Out port to charge external devices simultaneously with the batteries.

Easy travelling in focus




The travel adapter included in the set of **SERVOX® digital XL** guarantees easy and fast recharging of the **SERVOX® digital XL** in nearly every country in the world. The high quality materials and easy operation make it a reliable companion during your travels.

Individuality with sound frequency and quality




The **SERVOX® digital XL** speech aid offers the possibility to adjust and save a variety of individual settings. A more advanced fine-tuning of the sound parameters can be done with the help of our software, designed specifically for this purpose.

Variability of the sound nuances



The two screw caps of **SERVOX® digital XL** offer additional customization of the sound nuances. With the hard and soft screw cap you will find the right tone.

Variety in application



A diverse selection of mouth tube adapters gives an opportunity to use the **SERVOX® digital XL** speech aid immediately after a surgery or during a radiation therapy.

Mobile Companion



The **SERVOX® digital** belt pouch rounds off the **SERVOX® digital XL** set. It provides additional protection and is an ideal companion in everyday life. The integrated compartment can hold a second set of batteries as well as the additional battery holder as a backup.