HME STYLES & FEATURES GUIDE



	0.0	
I \	/ 1	
	r .	. 7













FEATURES	Blom-Singer® SpeakFree™ HME Hands- Free Valve with ClassicFlow® HME	Blom-Singer® SpeakFree™ HME Hands- Free Valve with EasyFlow® HME	Blom-Singer® ClassicFlow® HME	Blom-Singer® EasyFlow® HME	Blom-Singer® Day&Night ClassicFlow® HME	Blom-Singer® Day&Night EasyFlow® HME
Adjustable Hands Free Speech Valve	√	√				
Single-Use Device	√	√	✓	√	√	√
Easy Manual Occlusion for Voicing	✓	✓	✓	√	√	√
Breathability	GOOD	GOOD	BETTER	BETTER	BEST	BEST
Moisture Exchange	BETTER	GOOD	BETTER	GOOD	BETTER	GOOD
Hygroscopic Filter	√	√	√	√	√	√
Day and Night Wear	DAYTIME	DAYTIME	✓	✓	✓	✓
Compatible with all Blom-Singer Housings and StomaSoft Tube	√	√	1	√	√	√
Optional Safety Feature	CROSSBAR	CROSSBAR			TWIST TO LOCK	TWIST TO LOCK
MucusShield™ Protection			✓	√		
Made with Latex-free Materials	✓	✓	✓	√	√	✓
Non-reflective Surface Texture	✓	✓			√	√
Reference Number	BE 1090EZ 30 pk	BE 1090EF 30 pk	BE 1055EZ 30 pk	BE1055EF 30 pk	BE1091 30 pk	BE1092 30 pk

