

# Nasal Septal Perforation Prosthesis

**Blom-Singer®**

The voice of experience since 1978

## Two-piece magnetically coupled prosthesis

Magnetic  Closure 

A study of patients using a two-piece magnetically coupled prosthesis showed:

- No particular discomfort on insertion, removal, or replacement<sup>1</sup>
- Easier debris removal for ENTs<sup>1</sup>

### Choice of shape, size and thickness for a comfortable fit

Two medical-grade, silicone-covered neodymium magnets couple securely.

Multiple options – round or oval, 6 sizes each – help ensure an accurate fit:

- Less pistoning
- Less side-to-side movement

### Reduce or eliminate symptoms:

- Bleeding, crusting, whistling, epistaxis, infection<sup>2</sup>

Oval

Round

**INHEALTH®**  
TECHNOLOGIES

1300 005 279

**main**  
**MEDICAL**



# Nasal Septal Perforation Prosthesis

**Blom-Singer®**


The voice of experience since 1978

## OVAL



11 x 13 mm



13 x 17 mm



15 x 21 mm




Flange width is 5mm for all prostheses.



**Diagrams actual size**

## ROUND



5 mm



7 mm



10 mm



**INHEALTH®**  
TECHNOLOGIES

*We speak ENT*

1300 005 279

**main**  
**MEDICAL**