

INTRODUCING

INHEALTH
TECHNOLOGIES®

SpeakFree™ HME



New Blom-Singer® SpeakFree™ HME Hands Free Valve

It's simple. Your
HME and hands free
speaking valve, all
in one.



Blom-Singer® SpeakFree™ HME Hands Free Valve is the first heat and moisture exchange (HME) cartridge with a single-use fully integrated adjustable hands free valve.

Features:

- Choose your flow and go. Two simple choices for your breathability. ClassicFlow® or EasyFlow®.
- Truly innovative. One adjustable valve with hands free speech that adjusts from light to extra strong.
- Freedom of choice. Attaches to Blom-Singer® adhesive housings or a Blom-Singer® StomaSoft® Laryngectomy Tube.

TOTAL
Solution

SpeakFree™ HME




Blom-Singer® SpeakFree™ HME Hands Free Valve with ClassicFlow® HME

- Breathe healthy
- ClassicFlow® helps assure you breathe warmed, filtered and humidified air - every day - to reduce respiratory symptoms after a total laryngectomy. ClassicFlow® is designed for increased moisture retention and resistance.

Blom-Singer® SpeakFree™ HME Hands Free Valve with EasyFlow® HME

- Breathe easier at work or play
- With EasyFlow®, the foam filter is more porous than the ClassicFlow®. This allows you to breathe more freely to match your activity level and pulmonary needs.


One hands free valve is an alternative to multiple competitor devices.





One SpeakFree™ HME cartridge can be adjusted to do the work of multiple hands free competitor devices.

A SpeakFree™ HME all-in-one cartridge can be adjusted from as little as 45 liters of air to as much as 115 liters of air to attenuate the speaking valve.

Data on file.



Whatever your activity or lifestyle. A simple, disposable hand free HME that fits with your active life.

REFERENCE #	PRODUCT NAME	QTY
 BE 1090EF	Blom-Singer® SpeakFree™ HME Hands Free Valve with EasyFlow® HME <i>5 optional crossbars included</i>	30 per pack
 BE 1090EZ	Blom-Singer® SpeakFree™ HME Hands Free Valve with ClassicFlow® HME <i>5 optional crossbars included</i>	30 per pack

INHEALTH
TECHNOLOGIES®

Stay Connected with Us



©2021 InHealth Technologies®—Manufactured by Freudenberg Medical, LLC (210710.00)