


STERLINK Low Temperature PLASMA STERILIZER

Compatible Items

- Cranial pressure transducer cables
- Cryoprobes
- Defibrillator paddles
- Dopplers
- Electrocautery instruments
- Endoscopic instruments
- Esophageal dilators
- Fiberoptic lite cables
- Laryngoscope blades
- Laser handpieces, fibers, and accessories
- Metal instruments
- Ophthalmic lenses (diagnostic, magnifying)
- Patient lead cables
- Pigmentation handpieces
- Radiation therapy equipment
- Resectoscope / working elements and sheaths
- Rigid endoscopes
- Shaver handpieces
- Flexible endoscopes
- Stereotactic equipment and batteries
- Trocar sheaths
- Ultrasound probes
- Video cameras and couplers
- Micro instrument kit
- Laparoscopic Grasping Forceps
- Bipolar Forceps with cables
- Vessel Sealer
- Electrosurgical Pencils
- Endotracheal Tube
- Surgical Staples
- Micro Instrument Set
- Oxygen tubing
- Suction & Irrigation tube
- Oxygen mask
- Surgical tips
- Catheters
- Biopsy forceps




Compatible Materials

Metal	<ul style="list-style-type: none"> • Aluminum • Stainless Steel
Non-metal	<ul style="list-style-type: none"> • Acrylonitrile butadiene styrene (ABS) • Delrin • Ethyl vinyl acetate (EVA) • Fluorinated ethylene propylene (PFTE) • Glass, USP Type 1 borosilicate • High density polyethylene • Kraton • Latex • Polyethylene • Monel • Phenolic resin • Polycarbonate (PC) • Polyetherimide (PEI) • Polyethylene terephthalate (PET) • Polymethyl methacrylate (PMMA) • Polyphenyl sulphone • Polypropylene (PP) • Polystyrene • Polysulfone • Polytetrafluoroethylene (PTFE) • Polyurethane • Polyvinyl chloride (PVC) • Silicon

Non-Compatible Materials

Absorptive Materials (Cellulose)

- Wood
- Linen
- Paper
- Sponge



Liquids (Moisture)



Oil

Powders

