


vive® FOREARM CRUTCH

THANK YOU!

We sincerely appreciate your business. We strive to provide you with the best quality products at a great value.



COMFORTABLE GRIPS
& FOREARM REST



LIGHTWEIGHT



HEAVY DUTY



ADJUSTABLE HEIGHT

www.vivehealth.com


© Copyright and Distributed by Vive Health
8955 Fontana Del Sol Way, Naples, FL 34109
All Rights Reserved. Made in China.
Email: service@vivehealth.com



vivehealthusa


FIND CORRECT FOREARM CRUTCH HEIGHT

1. Let your arm hang naturally by your waist. Adjust the forearm crutch height so the handle touches your wrist.



CANE HEIGHT ADJUSTMENT

1. Press in the button snap.
2. Slide extension tube up or down to desired height. Button snap should secure cane into place.



WARNING

- Make sure the brass pin is locked into place before use.
- Weight Capacity: 250 LBS

GUARANTEE

You are protected by Vive Health's industry leading guarantees and customer service:



If you did not purchase through vivehealth.com, activate your guarantee at: www.vivehealth.com/register.

Product Code: MOB1014