

LEG COMPRESSION MACHINE

QUICK START GUIDE

1. Open the filter cover on the bottom of the unit and remove the screw.




2. Connect the power cord to a power outlet.

3. Securely connect the air hose to the pump.




4. Connect the air hose to each leg cuff, making sure to match the colored connectors with the correct color hose ending.



5. Insert your legs into the leg cuffs and zip both all the way up. Use the extenders if needed.



6. Press the Power button. Select your mode and time using the Mode and Time buttons. Press Start to start the pump.



For more information, please consult the User's Guide or visit us at ViveHealth.com

vive®

 For More Information, Visit:
<http://bit.ly/compressionpumpfaq>