

GOT MORE QUESTIONS?

Check out our list of Frequently Asked Questions at [vhealth.link/6dd78](https://vivehealth.link/6dd78) for helpful answers.



And if that doesn't answer your question, our customer service team would love to help! Feel free to connect with them by phone, e-mail, or chat on our website.



service@vivehealth.com



1-800-487-3808



vivehealth.com

vive®



TRANSFER BELT

Owner's Manual
RHB101E

vivehealth.com

OVERVIEW


We are constantly answering questions and recording helpful videos to make using your Vive Transfer Belt as easy as possible. Check out the included link and QR code to help you through the process.




To see all of the FAQ's in one place visit vhealth.link/6dd78

QUICK START GUIDE

1. Insert strap into the buckle through the end with the teeth.



2. Slide strap through the hole opposite the teeth.



QUICK RELEASE

To quickly release the strap, pull the release tab and slide the strap free.

