

THANK YOU!

We sincerely appreciate your business.
We strive to provide you with the best quality products
at a great value.


If you have any questions please contact us at

✉ service@vivehealth.com

☎ 1-800-487-3808
Monday - Friday 9am - 5pm EST

vive® WOODEN TRANSFER BOARD

User Manual
RHB1O37WOOD



www.vivehealth.com


© Copyright and Distributed by Vive Health
8955 Fontana Del Sol Way, Naples, FL 34109
All Rights Reserved. Made in China.
Email: service@vivehealth.com



vivehealthusa

HOW TO USE YOUR TRANSFER BOARD:

- 1) Place a portion of the board on the surface you are transferring from (1/3 suggested).
- 2) Place a portion of the board on the surface you are transferring to (1/3 suggested).
- 3) Support body weight on board and slowly move from one end to the other.



CAN BE USED ON:

- Wheel Chair
- Bed
- Chair

ACTIVATE YOUR GUARANTEE



You are protected by Vive Health's industry leading guarantees and customer service:

If you did not purchase through vivehealth.com, please register at vivehealth.com/register to validate your guarantee.

Product Code: RHB1O37WOOD