

Features & Benefits



1. The lightest power chair available on the market and ONLY available through Journey.
2. Weight is only 26 lbs without the battery
3. Powered by two 150w brushless motors, which work independently, enabling this agile power chair to turn on the spot.
4. 10Ah lithium battery, with extra battery storage place.
5. Joystick can easily be operated with one finger.
6. Joystick can function on left or right side of chair. Easy assembly to make the switch.
7. Electronic brakes instead of manual brakes make it safer when going up or down a slope.
8. Chair folds up in seconds and is perfect for travel, while also saving storage space.
9. USB charging port on the joystick for mobile phone.
10. Armrests easily lift up for easy transfer to and from the chair.

Specifications

SKU	08642	Battery Charge Time	6-8 hours
Overall Dimension	35.4" x 22.8" x 33.9"	Theoretical Travel Ranging	≥ 10 miles
Seat Height	18.1 in	Static Stability	≥ 6°
Seat Width	16.5 in	Dynamic Stability	≥ 3°
Seat Depth	15.7 in	Slope Performance	6°
Arm Height	9.8 in	Climbing Ability	≥ 3°
Seatback Height	15.7 in	Motor Type	24V/150W*2
Weight	26 lbs (without battery), 29.08 lbs (with battery)	Lithium Battery Type	24V/10Ah
Maximum Speed	≤ 2.8 mph	Max. Output Current of Controller	30A
Braking Distance	≤ 39.3 in	Max. Output Current of Charger	2A
Minimum Turning Radius	≤ 35.4 in	Front Wheel	Flat-Free Tire, O.D. 6.5 in
Weight Capacity	240 lbs	Rear Wheel	Flat-Free Tire, O.D. 10.6 in