

Features & Benefits



1. Among the lightest power chairs available on the market.
2. Powered by two 150W brushless motors which work independently, enabling this agile power chair to turn on the spot.
3. Electronic brakes instead of manual brakes make it safer when going up or down a slope.
4. Chair folds up in seconds and is perfect for travel while also saving storage space.
5. USB charging port on the joystick for mobile phone.
6. Armrests easily lift up for easy transfer to and from chair.
7. Battery pack design is easy to remove and charge.
8. Frame material is magnesium aluminum alloy.

Specifications

| | | | |
|------------------------|--|-----------------------------------|----------------------------|
| SKU | 08643 | Battery Charge Time | 6-8 hours |
| Overall Dimension | 36.22" x 22" x 35" | Theoretical Travel Ranging | ≥ 10 miles |
| Seat Height | 18 in | Static Stability | ≥ 6° |
| Seat Width | 16.5 in | Dynamic Stability | ≥ 3° |
| Seat Depth | 15.7 in | Slope Performance | 6° |
| Arm Height | 9.8 in | Climbing Ability | ≥ 3° |
| Seatback Height | 16.5 in | Motor Type | 24V/150W*2 |
| Weight | 34.83 lbs (without battery), 37.47 lbs (with battery) | Lithium Battery Type | 24V/6Ah |
| Maximum Speed | ≤ 2.8 mph | Max. Output Current of Controller | 35A |
| Braking Distance | ≤ 39.3 in | Max. Output Current of Charger | 2A |
| Minimum Turning Radius | ≤ 35.4 in | Front Wheel | Flat-Free Tire, O.D. 6 in |
| Weight Capacity | 240 lbs | Rear Wheel | Flat-Free Tire, O.D. 12 in |