


DeBug Beach Walker Standard Tall



Seat Dimensions
17w x 12.75d
Weight Capacity
350
Overall Weight
40

Perfect For People

5'6 to 6'4



Overall Dimensions High: 24w x 29.25d x 41.25h
Overall Dimensions Low (inches): 24w x 26.25d x 33h


Material Seat: PU+PE

Material Hand Grip: TPR

Wheel Size : 30cm



DeBug Beach Walker Petite Size



Seat Dimensions
14.5 w x 13.5 d
Weight Capacity
300
Overall Weight
40

Perfect For People
5'0 to 5'6



Overall Dimensions High (Inches): 24w x 26d x 34.54h

Overall Dimensions Low (Inches): 24w x 23.75d x 35h

Width 35"

Material Seat: PU+PE

Material Hand Grip: TPR

Wheel Size : 30cm

Shop Address
78 E Nine Mile Rd
Pensacola, FL 32534

office@debugmobiliyproducts.com

Corporate Address
Deming Designs, Inc.
1090 Cobblestone Dr,
Pensacola, FL 32514