

3 FUNCTION FULL & SEMI ELECTRIC HOME CARE BED



HALF RAIL



FULL RAIL



Hand Crank for HI/LO function on Semi Electric Bed



Mattress retainer holds mattress firmly in place

Sturdy 2 piece reinforced tubular deck

Built-in emergency hand crank lowers and raises the bed

4 x 3" castors - 2 locking and 2 swivel

Durable and washable bed ends

Bed Dimensions: 88" L x 36" W x 15" - 23 H

Side rails - Full or half rail optional

LOW POSITION: 15"

HIGH POSITION: 23"

HEAD ANGLE: 0 - 70 Degrees

FOOT ANGLE: 0 - 45 Degrees

WEIGHT CAPACITY: 450 Lbs.



MULTI FUNCTION PENDANT

- Raise and lower the head deck
- Raise and lower the foot deck
- Raise and lower the height of the bed - Full Electric

MODEL #: HCFE36 - FULL ELECTRIC

MODEL #: HCSE36 - SEMI ELECTRIC