



Echo

The Echo scooter has a light weight frame designed for the user on the go. The scooter disassembles into three pieces for portability and storage; with the heaviest part weighing just 35 lb. It has a height adjustable seat, allowing for improved comfort during driving operation. The Echo 3 also comes standard with a quick connect battery pack, solid tires and a basket to store personal items.

- + 250lb weight capacity
- + Quick connect battery pack
- + Removable battery for charging
- + Adjustable seat height

● Red ● Blue

Model Number: SL73

Technical Specifications

Overall Dimensions (L x W x H)	in	37" x 21" x 36"
Number, Size of Tires	Front, Rear	1, 8" (2, 8")
Suggested User Weight ¹	lb	Max. 250
Weight of Heaviest Part ²	lb	35
Weight ³	lb	70
Each Battery Module Weight	lb	15
Battery Size	Type x PCS	12V10AH x 2
Charger	Amp	Off board (2A)
Output Power of Motor	HP	0.45
Maximum Speed	MPH	3.8
Maximum Climbing Angle/	%(°)	10(6)
Safe Climbing Angle ⁴		
Range (Per Charge) ⁶	Mile	8
Turning Radius	in	32"
Ground Clearance	in	1.5"
Handle Bar	Type	T-Bar
Model #		SL73
Standard Seat	W x D x BH	SF171413B

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.