



## Operational warning for Glass Lined Tubing (GLT™)

### Important

#### **⚠ Warning**

Connections for GLT™ must not be made using metal compression couplings. For this reason, enclosed are 2 x OGF-8-8 ferrules.

### Information and support

Visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)

*Specifications are subject to change without notice.*

MN-0144-A\_RevC © Trajan Scientific Australia Pty Ltd 01/2017

[www.trajanscimed.com](http://www.trajanscimed.com)