




Instructions for syringe removal

Important


⚠ Caution

Be very careful not to break the foil to expose the next cavity with syringe until required.

Instructions



1. Push the backflange of the syringe (see figure below) to break the foil backing.



2. Once the foil is punctured at the backflange, the syringe can be removed by gently sliding it out. If required, the foil may be broken further by carefully running a thumb nail along the side of the cavity.

Information and support

Visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Specifications are subject to change without notice.

Analytical syringes manufactured by Trajan Scientific and Medical are intended for analytical and laboratory use only and are not intended or approved for use with food, including the production or packaging of food, nor medical or human in-vivo use.

MN-0334-S_RevD © Trajan Scientific Australia Pty Ltd 06/2020