

Instructions

Not all syringe part numbers are recognized by the PAL RTC instrument software. When using the SGE syringe part numbers below, enter the corresponding part number for instrument configuration.

SGE syringe part number	Description	Part number for instrument configuration
001871	5R-RTC/RSH-5.7/0.63C 5UL SYRINGE	001861
001873	5R-RTC/RSH-8.5/0.63C 5UL SYRINGE	001863
001875	5R-RTC/RSH-5.7/0.47C 5UL SYRINGE	001865
001877	5R-RTC/RSH-8.5/0.47C 5UL SYRINGE	001867
002790	10R-RTC/RSH-GT-MS1 10UL DIAMOND MS SYRINGE	002862
002852	10F-RTC/RSH-GT-5.7/0.63(0.24)C 10UL SYRINGE	002862
002854	10F-RTC/RSH-GT-8.5/0.72LC 10UL SYRINGE	002864
002869	10F-RTC/RSH-8.5/0.47 BV 10UL SYRINGE	002867
002871	10R-RTC/RSH-5.7/0.63C 10UL SYRINGE	002861
002872	10R-RTC/RSH-GT-5.7/0.63C 10UL SYRINGE	002862
002873	10R-RTC/RSH-8.5/0.63C 10UL SYRINGE	002863
002874	10R-RTC/RSH-GT-8.5/0.63C 10UL SYRINGE	002864
002875	10R-RTC/RSH-5.7/0.47C 10UL SYRINGE	002865
002876	10R-RTC/RSH-GT-5.7/0.47C 10UL SYRINGE	002866
002877	10R-RTC/RSH-8.5/0.47C 10UL SYRINGE	002867
002878	10R-RTC/RSH-GT-8.5/0.47C 10UL SYRINGE	002868
002880	10R-RTC/RSH-GT-5.7/0.72LC 10UL SYRINGE	002862
003790	25R-RTC/RSH-GT-MS1 25UL DIAMOND MS SYRINGE	003862
003852	25F-RTC/RSH-GT-5.7/0.63(0.24)C 25UL SYRINGE	003862
003854	25F-RTC/RSH-GT-8.5/0.72LC 25UL SYRINGE	003864
004790	50R-RTC/RSH-GT-MS1 50UL DIAMOND MS SYRINGE	004862
005790	100R-RTC/RSH-GT-MS1 100UL DIAMOND MS SYRINGE	005862
005852	100F-RTC/RSH-GT-5.7/0.63(0.24)C 100UL SYRINGE	005862
005854	100F-RTC/RSH-GT-8.5/0.72LC 100UL SYRINGE	005864
006854	250F-RTC/RSH-GT-8.5/0.72LC 250UL SYRINGE	006864
008151	1MR-RTC/RSH-GT-5.7/0.72LC 1ML SYRINGE	008150
008951	10MR-RTC/RSH-GT-5.7/1.07LC 10ML SYRINGE	008950

Information and support

Visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Specifications are subject to change without notice.