

## Separations | HPLC guard column ProteCol®

Turn imagination into results

ProteCol® guard columns are recommended to protect the analytical column, and ensure it performs consistently.

### Features and benefits

The unique design of the guard column couples directly to the analytical column. The guard column is designed to fit into the back of a PEEK fingertight fitting (provided with the guard column). No further unions are required.

### Recommended applications

We recommend the use of ProteCol Guard columns to protect the analytical column, and ensure it performs consistently.

### Product specifications

Available in C18 and C8 phases.

For more information about this product visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)

