

## Connections | Graphite/Vespel® ferrules

Part number	Part description and detail
072694	15% Graphite/85% Vespel mini union ferrule with no hole, PK10
072695	15% Graphite/85% Vespel mini union ferrule with 0.3 mm hole, PK10
072696	15% Graphite/85% Vespel mini union ferrule with 0.4 mm hole, PK10
072697	15% Graphite/85% Vespel mini union ferrule with 0.5 mm hole, PK10
072698	15% Graphite/85% Vespel mini union ferrule with 0.8 mm hole, PK10
072671	15% Graphite/85% Vespel ferrule with 0.8 mm hole for 1/8" fitting, PK10
072674	15% Graphite/85% Vespel ferrule with 2 x 0.8 mm holes for 1/8" fitting, PK10
072673	15% Graphite/85% Vespel ferrule with 1/16" hole for 1/4" fitting, PK10
<b>For use on Agilent GC injectors and detectors at atmospheric pressure e.g. FID</b>	
073105	15% Graphite/85% Vespel short ferrule with no hole, PK10
073107	15% Graphite/85% Vespel short ferrule with 0.3 mm hole, PK10
073109	15% Graphite/85% Vespel short ferrule with 0.4 mm hole, PK10
073111	15% Graphite/85% Vespel short ferrule with 0.5 mm hole, PK10
073113	15% Graphite/85% Vespel short ferrule with 0.8 mm hole, PK10
<b>For use on Agilent GCMS interface</b>	
072663	15% Graphite/85% Vespel ferrule with 0.4 mm hole, PK10
072654	15% Graphite/85% Vespel ferrule with 0.5 mm hole, PK10
072655	15% Graphite/85% Vespel ferrule with 0.8 mm hole, PK10
<b>For use on Bruker/Varian GC injectors and detectors at atmospheric pressure e.g. FID</b>	
072663	15% Graphite/85% Vespel ferrule with 0.4 mm hole, PK10
072654	15% Graphite/85% Vespel ferrule with 0.5 mm hole, PK10
072655	15% Graphite/85% Vespel ferrule with 0.8 mm hole, PK10
<b>For use on PerkinElmer GC injectors and detectors at atmospheric pressure e.g. FID with 1/16" nut</b>	
072663	15% Graphite/85% Vespel ferrule with 0.4 mm hole for 1/16" fitting, PK10
072654	15% Graphite/85% Vespel ferrule with 0.5 mm hole for 1/16" fitting, PK10
072655	15% Graphite/85% Vespel ferrule with 0.8 mm hole for 1/16" fitting, PK10
<b>For use on PerkinElmer GC injectors and detectors at atmospheric pressure e.g. FID with 1/8" nut</b>	
0726703	15% Graphite/85% Vespel ferrule with 0.4 mm hole for 1/8" fitting, PK10
0726702	15% Graphite/85% Vespel ferrule with 0.5 mm hole for 1/8" fitting, PK10
072671	15% Graphite/85% Vespel ferrule with 0.8 mm hole for 1/8" fitting, PK10
<b>For use on Shimadzu QP5000-I standard MS interface</b>	
0726563	15% Graphite/85% Vespel ferrule with 0.4 mm hole, PK10
0726564	15% Graphite/85% Vespel ferrule with 0.5 mm hole, PK10
<b>For use on Shimadzu QP5000-II/5050 standard MS interface</b>	
0726561	15% Graphite/85% Vespel ferrule with 0.4 mm hole, PK10
0726562	15% Graphite/85% Vespel ferrule with 0.5 mm hole, PK10
0726547	15% Graphite/85% Vespel ferrule with 0.8 mm hole, PK10
<b>For use on Shimadzu QP5000/5050 wide bore MS interface, QP2010 injector and QP2010 standard MS interface</b>	
072663	15% Graphite/85% Vespel ferrule with 0.4 mm hole, PK10
072654	15% Graphite/85% Vespel ferrule with 0.5 mm hole, PK10
072655	15% Graphite/85% Vespel ferrule with 0.8 mm hole, PK10

## Connections | Graphite/Vespel® ferrules continued

Part number	Part description and detail
<b>For use on Thermo Scientific TRACE 1300 Series GC at atmospheric pressure e.g. FID</b>	
073105	15% Graphite/85% Vespel short ferrule with no hole, PK10
073107	15% Graphite/85% Vespel short ferrule with 0.3 mm hole, PK10
073109	15% Graphite/85% Vespel short ferrule with 0.4 mm hole, PK10
073111	15% Graphite/85% Vespel short ferrule with 0.5 mm hole, PK10
073113	15% Graphite/85% Vespel short ferrule with 0.8 mm hole, PK10
<b>For use on Thermo Scientific GCMS interface</b>	
072663	15% Graphite/85% Vespel ferrule with 0.4 mm hole, PK10
072654	15% Graphite/85% Vespel ferrule with 0.5 mm hole, PK10
072655	15% Graphite/85% Vespel ferrule with 0.8 mm hole, PK10
<b>Standard ferrule solutions - 1/16" fitting specification</b>	
0726632	15% Graphite/85% Vespel ferrule with 0.2 mm hole for 1/16" fitting, PK10
0726591	15% Graphite/85% Vespel Vg/Micromass ferrule with 0.3/0.4 mm hole for 1/16" fitting, PK10
0726650	15% Graphite/85% Vespel ferrule with 3 x 0.3 mm holes for 1/16" fitting, PK10
0726641	15% Graphite/85% Vespel ferrule with 1 x 0.4 mm hole and 1 x 0.3 mm hole for 1/16" fitting, PK10
0726640	15% Graphite/85% Vespel ferrule with 1 x 0.5 mm hole and 1 x 0.3 mm hole for 1/16" fitting, PK10
0726642	15% Graphite/85% Vespel ferrule with 1 x 0.5 mm hole and 1 x 0.4 mm hole for 1/16" fitting, PK10
<b>Standard ferrule solutions - 1/8" fitting specification</b>	
0726703	15% Graphite/85% Vespel ferrule with 0.4 mm hole for 1/8" fitting, PK10
0726702	15% Graphite/85% Vespel ferrule with 0.5 mm hole for 1/8" fitting, PK10
072674	15% Graphite/85% Vespel ferrule with 2 x 0.8 mm holes for 1/8" fitting, PK10
0726701	15% Graphite/85% Vespel ferrule with 1/16" hole for 1/8" fitting, PK10
0726700	15% Graphite/85% Vespel ferrule with 2.4 mm hole for 1/8" fitting, PK10
<b>Standard ferrule solutions - 1/4" fitting specification</b>	
0726660	15% Graphite/85% Vespel ferrule with 0.5 mm hole for 1/4" fitting, PK10
072673	15% Graphite/85% Vespel ferrule with 1/16" hole for 1/4" fitting, PK10
0726673	15% Graphite/85% Vespel ferrule with 6 mm hole for 1/4" fitting, PK10

Vespel® is registered by DuPont.

For more information about this product visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)