



Separations | ProDx columns

Optimized for HDX workflows

Part number	Part description and detail
Protease	
PDX.PP01-F32	30 mm x 2.1 mm ID ProDx Protease (Pepsin) bio-inert stainless steel glass lined column
Trap	
PDX.T800	20 mm x 530 µm ID ProDx C8 Analog* trap column PK3 (without fittings)
PDX.T8L02-32	20 mm x 530 µm ID ProDx C8 Analog* trap column PK3 (with 1/16" to 10-32 UNF thread (long) ProteCol® PEEK Fingertight HPLC fittings PK6)
PDX.T8L02-40	20 mm x 530 µm ID ProDx C8 Analog* trap column PK3 (with 1/32" to 6-40 UNF thread ProteCol PEEK Fingertight HPLC fittings PK6)
PDX.T8H02-32	20 mm x 530 µm ID ProDx C8 Analog* trap column PK3 (with 1/16" to 10-32 UNF thread PEEKLok™ UHPLC fittings (6 x adaptor ferrules, 2 x hex head fittings))
PDX.T8H02-40	20 mm x 530 µm ID ProDx C8 Analog* trap column PK3 (with 1/32" to 6-40 UNF thread PEEKLok UHPLC fittings (6 x ferrules, 2 x hex head fittings))

**This product contains POROS® material used under license for Sample Prep purposes from Life Technologies Corporation, a subsidiary of Thermo Fisher Scientific, Inc. The POROS® Trademark is owned by Thermo Fisher Scientific, Inc.*

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com