

Training cards

Lesson 1

- 1. Creating a write cell
- 2. Changing a picture
- 3. Editing a picture
- 4. Moving a cell
- 5. Editing the grid layout
- 6. Changing cell appearance
- 7. Creating jump cells
- 8. Chat and writing cells

Lesson 2

- 1. Introducing commands
- 2. Word lists
- 3. Adding grids
- 4. Cell accessibility
- 5. Adding sound to a cell
- 6. Setting up your accounts
- 7. Send to Grid Player
- 8. Remote editing

Lesson 3

- 1. Settings Grid set
- 2. Settings Email accounts
- 3. Settings Speech
- 4. Settings Mouse or head pointer
- 5. Settings Touch
- 6. Settings Eye gaze
- 7. Settings Switch
- 8. Settings Mobile and cell phone



Visit our website at www.thinksmartbox.com to access training videos, guides and support.





Edit mode keyboard shortcuts

F11	-	Start/Finish editing
Ctrl + X Ctrl + C Ctrl + V Ctrl + Z Ctrl + S Ctrl + A Ctrl + D Ctrl + D Ctrl + P Ctrl + J Ctrl + F Ctrl + W	- - - - - - - -	Cut Copy Paste Undo Save Select all Deselect all Change cell picture Create jump Follow jump Create write cell
Arrows Enter Insert Delete	- - -	Navigate cells (+ shift to select) Create cell Add command (or create cell if blank) Clear cell
F1 F2 F3 F4 F5 F10 F12	- - - - -	Help Change cell label Find picture Edit picture Remove picture Settings Show/Hide menu bar in Grid Viewer
Alt	-	Bring up Keytips to navigate the ribbon menus with keyboard letters

