# Included with your Grid Pad<sup>™</sup> 13:



## Hardware, Buttons, & Features

#### Power

Press this button to turn on and off the device.

## ·Þ

#### **USB** Ports

Use ports for connecting accessibility options such as a mouse, keyboard, or switch interface.

### • Headphone/Microphone


Audio input/output for sound feedback, voice and message banking and more.

🜒 Volume

Press to adjust volume levels.

## Charging the Grid Pad 13

Charge device at night or when not in use. Please only use TTMT-approved chargers with your device.



Grid Pad 13





