



Certificate of Quality

TipOne® Filter Tips

Article No.	10 µL	10/20 µL XL	200 µL	300 µL	1000 µL	1000 µL XL
Filter Tips	9104020341	9104020357	9104020365	9104020497	9104020872	9104020732

Eppendorf certifies that this product is produced in a fully controlled environment that ensures the product is free from detectable RNase, DNase, DNA and Endotoxins. Selected samples of product from this production environment have been tested for the following parameters on a recurring basis.

- **RNase testing method:** RNA Digestion
Limits: $< 1 \times 10^{-9}$ Kunitz Units
- **DNase testing method:** DNA Digestion
Limits: $< 1 \times 10^{-6}$ Kunitz Units
- **Human DNA testing method:** Polymerase Chain Reaction
Limits: $< 2 \times 10^{-9}$ mg
- **Endotoxin testing method:** Kinetic-chromogenic LAL-Test
Limits: < 0.05 EU/mL

TipOne® Filter Tips have been treated by e-beam irradiation and the process validated according to ISO 11137 'Sterilization of health care products'. The validation was performed by using 15 kGy as the minimum dose with SAL (sterility assurance level) of 10^{-6} .

Hamburg, October 2023

Bhavna Mungur
Head of Quality

Thomas Hengstmann
Head of Global Quality Operations

Eppendorf SE, Barkhausenweg 1, 22339 Hamburg, Germany
Phone: + 49 40 53 801-0

www.excella.bio