



## Certificate of Quality

### TipOne® Pipette Tips

Article No.	10 µL	10/20 µL XL	200 µL	300 µL	1000 µL	1250 µL XL
<b>Rack</b>	9104020317	9104020279	9104020406	9104020470	9104020578	9104020716
<b>Refill</b>	9104020333	9104020295	9104020422	9104020462	9104020551	9104020740
<b>Bulk</b>	9104020309	9104020260	9104020449	9104020454	9104020543	9104020708

Eppendorf certifies that these products are produced in a fully controlled environment that ensures the products are free from detectable RNase, DNase, DNA and Endotoxins. Selected samples of product from this production environment have been tested for the following parameters on a recurring basis.

- **RNase testing method:** RNA Digestion  
**Limits:**  $< 1 \times 10^{-9}$  Kunitz Units
- **DNase testing method:** DNA Digestion  
**Limits:**  $< 1 \times 10^{-6}$  Kunitz Units
- **Human DNA testing method:** Polymerase Chain Reaction  
**Limits:**  $< 2 \times 10^{-9}$  mg
- **Endotoxin testing method:** Kinetic-chromogenic LAL-Test  
**Limits:**  $< 0.05$  EU/mL

Hamburg, November 2023

Bhavna Mungur  
Head of Quality

Thomas Hengstmann  
Head of Global Quality Operations

Eppendorf SE, Barkhausenweg 1, 22339 Hamburg, Germany  
Phone: + 49 40 53 801-0

[www.excella.bio](http://www.excella.bio)