



Manufacturing Chemists Since 1953



*For Qualitative
Determination of
Fecal Occult Blood*



www.ecamco.com

Made and Printed in the U.S.A.

GUAIAIC-TABS®



For *In Vitro* Diagnostic Use

- *For Qualitative Determination of Fecal Occult Blood* •
- *Quick, Simple Procedure* • *Results in 2 minutes* •
- *100 Filter Paper Squares Provided* •

1003, Guaiac-Tabs are made of refined and purified Guaiac Gum and are used to determine the presence of occult blood in a fecal sample. Each box contains 100 tablets and 100 filter paper squares for your convenience.

This simple, reliable test requires only a tablet, a filter paper square, glacial acetic acid and 3% hydrogen peroxide to complete the procedure in only 2 minutes. Development of a blue color on the filter paper may indicate the presence of blood.

Available in the following size

1003	100 Tablets	100 Tablets with 100 Filter Paper Squares
------	-------------	-------------------------------------------

Available from your laboratory supply distributor or directly from:

Cambridge Diagnostic Products, Inc.
6880 NW 17th Avenue
Fort Lauderdale, FL 33309

For additional information:
techinfo@ecamco.com
954.971.4040 or 800.525.6262
www.ecamco.com



Manufacturing Chemists Since 1953



www.ecamco.com

Made and Printed in the U.S.A.