







Manufacturing Chemists Since 1953

# For Automated and Manual Staining of Blood Smears



[www.ecamco.com](http://www.ecamco.com)

Made and Printed in the U.S.A.



# CAMCO<sup>®</sup>

## Wright Stain, Fixative and Buffer



For In Vitro Diagnostic Use

- *For Automated and Manual Staining of Blood Smears* •
- *Classic, Traditional Textbook Coloration* •

The unique combination of these products has been developed for the use in automated and manual staining of blood smears for hematological studies. These special 5-step processes (featured on the Camco<sup>®</sup> Wright Stain label) provides stained blood smears of classic tradition with textbook colors utilizing our Camco<sup>®</sup> Wright Stain, Camco<sup>®</sup> Fixative and Camco<sup>®</sup> Wright Stain Buffer pH 7.15 or pH 6.8

**400, CAMCO<sup>®</sup> Wright Stain** is highly refined, modified and purified alcoholic Wright Stain solution. This specially formulated stain makes analysis of the smear much easier.

**401, CAMCO<sup>®</sup> Fixative** is methanol-based fixative to set the blood smears as the first step in the 5-step processes. **Also may be used for any fixing purposes where Methanol fixing is needed.**

**402, CAMCO<sup>®</sup> Wright Stain Buffer pH 7.15.** Each 4.4 g vial is enough to prepare 1 gallon of pH 7.15 buffer solution for use in the 5-step staining process. The buffer is helpful with the development of the correct, clear nuclear color scheme.

**412, CAMCO<sup>®</sup> Wright Stain Buffer pH 6.8** Each 4.4 g vial is enough to prepare 1 gallon of pH 6.8 buffer solution for use in the 5-step staining process. The buffer is helpful with the development of the correct, clear nuclear color scheme.

Available in the following sizes

400	1 gallon / 3.785 L	Camco Wright Stain
401	1 gallon / 3.785 L	Fixative
402	Ten - 4.4 g	Vials (Box) 7.15 pH buffer
412	Ten - 4.4 g	Vials (Box) 6.8 pH buffer

Available from your laboratory supply distributor or directly from:

**Cambridge Diagnostic Products, Inc.**  
6880 NW 17<sup>th</sup> Avenue  
Fort Lauderdale, FL 33309

For additional information:  
[techinfo@ecamco.com](mailto:techinfo@ecamco.com)  
954.971.4040 or 800.525.6262  
[www.ecamco.com](http://www.ecamco.com)



Manufacturing Chemists Since 1953



[www.ecamco.com](http://www.ecamco.com)

Made and Printed in the U.S.A.