



Manufacturing Chemists Since 1953



A Three-Step Significantly Superior Staining System



www.ecamco.com

Made and Printed in the U.S.A.

CAMCO[®] STAIN PAK



For *In Vitro* diagnostic use

- *Approximating results obtained from either Wright or Wright-Giemsa Stain* •
- *Perfectly differentiated color scheme* • *Easy to read results, every time* •
- *Quick* • *Simple* • *Bright* • *Clean* •

702-A, Fixative Solution is a methanol based fixative with added dye stuff.

702-B, Solution I and **702-C, Solution II** are two water-based solutions, separating the polychrome dyestuff into its base components. Staining can easily be adjusted to increase coloration to more red or blue. In addition to the highly purified dyestuff and buffers, **Camco Stain Pak** contains important electrolytes to enhance the staining of the nuclei and its constituents, offering you an even better analytical view.

Available in the following sizes

702	One Pak includes the following: three bottles each, 16 fl. oz. / 473 mL	Camco Stain Pak Fixative Solution I Solution II
702-A	1 gallon / 3.785 L	Fixative
702-B	1 gallon / 3.785 L	Solution I Red
702-C	1 gallon / 3.785 L	Solution II Blue

Available from your laboratory
supply distributor or directly from:

Cambridge Diagnostic Products, Inc.
6880 NW 17th Avenue
Fort Lauderdale, FL 33309

For additional information:
techinfo@ecamco.com
954.971.4040 or 800.525.6262
www.ecamco.com



Manufacturing Chemists Since 1953



www.ecamco.com

Made and Printed in the U.S.A.