

# Glaucoma Punches, Rongeurs


## Glaukom- und Knochenstanzen

Rongeurs 65  
Glaucoma Punches 66



G-18220


BEYER  
**Bone Nibbling Rongeur**  
2 mm x 3 mm bite  
overall length 140 mm




Scale  
2:5

G-18225



KERRISON  
**Bone Nibbling Rongeur**  
5 mm x 3.5 mm bite  
overall length 160 mm




Scale  
2:5



☒-18238 KELLY  
**Descemet Membrane Punch**  
for Glaucoma  
19 Gauge/ 1.0 mm  
0.9 mm x 0.45 mm bite  
overall length 122 mm



☒-18239 LUNTZ - DODICK  
**Trabeculectomy Punch**  
for Glaucoma  
with 360° rotatable cutting head  
1 mm x 0.5 mm bite  
overall length 115 mm



☒-18235 GREHN  
**Micro Trabeculectomy Punch**  
for Glaucoma  
20 Gauge/ 0.9 mm  
0.7 mm x 0.3 mm bite  
overall length 90 mm

