

Geuder® **09**
Forceps

Pinzetten

- Cilia, Foreign Body 87
- Conjunctiva 88
- Suturing 89
- Keratoplasty 91
- Suturing, Tying 91
- Fragment, Jewelers 94
- Fixation, Dressing, Tissue 95
- Iris 95
- Capsulorhexis 97
- DMEK, DSAEK, LASIK 101
- Lenticule 102

- Nucleus Prechoppers,
Implantation 103
- IOL 104
- Superior Rectus 107
- Scleral Plug, Trocar 109
- Vitreoretinal 110
- Chalazion, Ptosis 112
- Mosquito 113
- Towel, Serrefine 114



Geuder® **09 Forceps**
Cilia, Foreign Body

G-18720

HENRY
Cilia Forceps
straight jaws
smooth, blunt, squared tips
overall length 100 mm



G-18730

GRADLE
Cilia Forceps
curved jaws
smooth, blunt, squared tips
overall length 98 mm



G-18750

BARRAQUER
Cilia Forceps
lightly angled
65 mm platform
overall length 109 mm



G-18740

BEER
Cilia Forceps
3 mm platform
rounded tips 45 mm wide
overall length 78 mm



G-18850

SHAAF
Cilia and Foreign Body Forceps
grooved tips
curved and pointed
finger grip handle
overall length 96 mm



G-31620

BARRAQUER
Cilia and Suture Forceps
grooved tips, 7 mm platform
curved and pointed
finger grip handle
overall length 112 mm




Geuder® 09 Forceps
Cilia, Foreign Body

GS-32740 HAMBURG MODEL
Cilia and Foreign Body Forceps
tapered fine tips
overall length 90 mm




Geuder® 09 Forceps
Conjunctiva


GS-31230 KREFELD MODEL
Conjunctiva Forceps
standard
1 x 2 teeth
overall length 100 mm



GS-31235 FREIBURG MODEL
Conjunctiva Forceps
delicate
1 x 2 teeth
overall length 90 mm




GS-18810 DRAEGER
Conjunctiva Forceps
12 mm wide
2 atraumatic fixation knobs
overall length 97 mm




Geuder® 09 Forceps
Suturing


CASTROVIEJO
Suturing Forceps
1 x 2 teeth
6 mm tying platform
overall length 100 mm
GS-19001 0.28 mm, teeth 0.12 mm
GS-19005 0.3 mm
GS-19010 0.6 mm
GS-19015 0.9 mm
GS-19020 1.2 mm



CASTROVIEJO
Suturing Forceps
new model
1 x 2 teeth
6 mm tying platform
overall length 100 mm
GS-19022 0.28 mm, teeth 0.12 mm
GS-19023 0.3 mm
GS-19024 0.5 mm
GS-19025 0.9 mm
GS-19026 1.2 mm




GS-19008 MUNICH MODEL
Suturing Forceps
1 x 2 teeth, 0.35 mm
55 mm tapered tying platform
overall length 100 mm



Geuder® 09 Forceps
Suturing

ZUERICH MODEL
Suturing Forceps
delicate
1 x 2 teeth, 0.12 mm
5 mm tying platform
overall length 75 mm
straight
angled


☒-18965
☒-18966



Geuder® 09 Forceps
Colibri Style


☒-18950

Colibri Forceps
1 x 2 teeth, 0.12 mm
5 mm tying platform
overall length 75 mm




☒-18960

Micro Colibri Forceps
delicate
1 x 2 teeth, 0.12 mm
5 mm tying platform
overall length 77 mm



☒-15134

Micro Colibri Forceps
extra delicate
1 x 2 teeth, 0.12 mm
4 mm tying platform
overall length 76 mm



☒-18975


BARRAQUER
Suture Removing Forceps
very fine pointed tips
4 mm tying platform
overall length 77 mm



Geuder® 09 Forceps
Keratoplasty


☒-18972

THOMAS NEUHANN
Keratoplasty Forceps
for corneal suturing and fixation
two prongs with 1 x 2 teeth
0.8 mm prong spread
overall length 80 mm




☒-18970

REMKY
Keratoplasty Forceps
two prongs with 1 x 1 teeth, angled
1.7 mm prong spread
for corneal suturing and fixation
overall length 80 mm



☒-38623


REMKY
DMEK Keratoplasty Forceps
modified by SZURMAN
two prongs with 1 x 1 teeth, straight
1.7 mm prong spread
33 mm shafts
to hold graft during preparation
overall length 89 mm



Geuder® 09 Forceps
Suturing, Tying


☒-18940
☒-18945

PAUFIQUE
Universal Suturing Forceps
1 x 2 teeth
6 mm tying platform
overall length 88 mm
standard
delicate



☒-19045
☒-19046

PAUFIQUE
Tying Forceps
with platform
overall length 88 mm
standard, 6 mm platform
delicate, 5 mm platform




Geuder® 09 Forceps
Suturing, Tying


GS-19040
TUEBINGEN MODEL
Tying Forceps
curved, smooth jaws
6 mm tying platform
overall length 99 mm




GS-19035
TUEBINGEN MODEL
Tying Forceps
delicate, smooth jaws
4 mm tying platform
overall length 84 mm




GS-19034
BONN MODEL
Tying Forceps
extra delicate, smooth jaws
3 mm tying platform
overall length 73 mm



GS-19032
Tying Forceps
angled, 7.5 mm
6 mm tying platform
overall length 99 mm




GS-31625
DARDENNE
Tying Forceps
angled, 6.5 mm
6 mm tying platform
round handle with guide pin
overall length 102 mm




Geuder® 09 Forceps
Suturing, Tying

GS-19038
BANGERTER
Tying Forceps
9.5 mm platform
longitudinal V-shape groove
overall length 100 mm




GS-19048
Tying Forceps
utility serrated jaws
1 mm wide
overall length 106 mm



GS-32947
HATTENBACH
Hybrid Tying Forceps
angled, for the anterior chamber
intraocular 1.5 mm micro tying tips
25 Gauge/0.5 mm
overall length 98 mm

for use in combination with
HATTENBACH Hybrid Needle Holder
GS-32949



Geuder® 09 Forceps
Fragment, Jewelers

GS-31623
KELMAN - McPHERSON
Fragment Forceps
angled, smooth jaws
7.5 mm platform
overall length 100 mm



GS-33941
KELMAN - McPHERSON
Fragment Forceps
angled, fine grooved jaws,
for Femto-LASIK to retrieve lenticule
during SMILE procedures
7.5 mm platform
overall length 100 mm



GS-31622
GS-31624
GASKIN
Fragment Forceps
angled
with tying platform
9 mm platform, overall length 107 mm
11 mm platform, overall length 100 mm



GS-37700
Jewelers Type Forceps
standard, straight
overall length 120 mm



GS-37705
Jewelers Type Forceps
delicate, straight
overall length 115 mm



GS-37710
Jewelers Type Forceps
extra delicate, straight
overall length 115 mm




Geuder® 09 Forceps
Fixation, Dressing, Tissue


GS-19050
GRAEFE
Fixation Forceps
without lock
4.5 mm wide
8 interlocking teeth
overall length 110 mm



GS-18780
GS-18781
GS-18782
GS-18783
GS-18784
Dressing Forceps
serrated jaws
overall length 100 mm




GS-18790
GS-18791
GS-18792
GS-18793
GS-18794
Tissue Forceps
1 x 2 teeth
overall length 100 mm




Geuder® 09 Forceps
Iris

GS-18855
GS-18856
Iris Forceps
serrated, delicate
overall length 73 mm
straight
curved



GS-18908
GS-18911
BONN MODEL
Iris Forceps
serrated, extra delicate
overall length 73 mm
straight
angled at 8 mm



Geuder® 09 Forceps


Iris

BONN MODEL

Iris Forceps

1 x 2 teeth, 0.12 mm
extra delicate
overall length 73 mm
straight
angled at 6 mm
angled at 8 mm

- ☒ -18900
- ☒ -18905
- ☒ -18910




BONN MODEL

Iris Forceps

long, straight
1 x 2 teeth, 0.12 mm
overall length 102 mm
extra delicate
delicate, with 45 mm tying platform

- ☒ -18845
- ☒ -18890




SWISS MODEL

Iris Forceps

curved, 1 x 2 teeth
overall length 100 mm


- ☒ -18840



Iris Forceps

straight, delicate
1 x 2 teeth
overall length 72 mm



- ☒ -18857



Iris Forceps

straight, strong
1 x 2 teeth
overall length 72 mm

- ☒ -31953




SWISS MODEL

Iris Forceps

straight, extra delicate
1 x 2 teeth
overall length 70 mm

- ☒ -32240



Geuder® 09 Forceps


Iris

- ☒ -18897

ELZE

Iris Forceps

extra delicate
1 x 2 teeth
overall length 78 mm




- ☒ -18913

THOMAS NEUHANN

Iris Forceps

curved shaft
1 x 2 teeth
23 Gauge/0.6 mm
overall length 95 mm



Geuder® 09 Forceps

Capsulorhexis

- ☒ -31307

UTRATA

Capsulorhexis Forceps

strong, with heel stop to
protect the iris
interlocking micro serrated tips
12 mm straight shafts, angled
overall length 103 mm



- ☒ -31308

Capsulorhexis Forceps

round handle with guide pin
extra delicate with
interlocking micro serrated tips
13 mm curved shafts, angled
overall length 105 mm



- ☒ -31300

MASKET


Capsulorhexis Forceps

delicate, with heel stop to
protect the iris
interlocking micro serrated tips
12 mm curved shafts, angled
overall length 103 mm




Geuder® 09 Forceps
Capsulorhexis


☒-31304 **Capsulorhexis Forceps**
new model, with guide pin
extra delicate with
interlocking micro serrated tips
13 mm curved shafts, angled
overall length 103 mm




☒-31301 **Capsulorhexis Forceps**
extra delicate
with smooth triangular grasping tips
13 mm curved shafts, angled
overall length 98 mm




☒-31305 **Capsulorhexis Forceps**
extra delicate with
interlocking micro serrated tips
13 mm curved shafts, angled
overall length 98 mm



☒-31303 **H.-R. KOCH
Capsulorhexis Forceps**
extra delicate with
interlocking micro serrated tips
13 mm curved shafts
1 mm marks from tip to bent
overall length 98 mm




☒-31299 **Capsulorhexis Forceps**
short handle, extra delicate
cross serrated grasping tips
12 mm straight shafts, angled
overall length 83 mm




Geuder® 09 Forceps
Capsulorhexis


☒-32926 **AUFFARTH
MICS Double Cross Action
Capsulorhexis Forceps**
interlocking micro serrated tips
12 mm curved shafts, angled 130°
with overlap system
prevents inadvertent iris capture
for 1.6 mm incisions
overall length 100 mm




☒-32928 **INAMURA
MICS Capsulorhexis Forceps**
cross action
interlocking micro serrated tips
105 mm curved shafts tip to pivot
angled 130°
overall length 103 mm



☒-32935 **H.-R. KOCH
MICS Capsule Forceps**
short handle
22 Gauge/0.7 mm curved shaft with
micro grasping tips
overall length 98 mm




☒-32938 **H.-R. KOCH
MICS Capsule Forceps**
22 Gauge/0.7 mm curved shaft with
micro grasping tips
6 mm titanium handle, blue-colored
overall length 144 mm




Geuder® 09 Forceps
Micro Incision



▣-32936 TOBIAS NEUHANN
Micro Forceps
multifunction asymmetrical tips with grasping platform
23 Gauge/ 0.6 mm straight shaft
6 mm titanium handle, blue-colored
overall length 155 mm



IKEDA
Capsulorhexis Forceps
interchangeable tips
23 Gauge/ 0.6 mm curved shaft
overall length 40 mm
▣-38244 textured, micrograsping tips "BEAK"
▣-38242 textured, micrograsping tips
for use with JAPAN MODEL Handle
▣-38240




▣-38244 "BEAK" ▣-38242




▣-38240 JAPAN MODEL
Handle
for capsule scissors and forceps
interchangeable tips
6 mm, titanium
overall length 90 mm



▣-32944 HATTENBACH
Capsulorhexis Forceps
25 Gauge/ 0.5 mm angled shaft with micro grasping tips
overall length 88 mm




▣-38624 SZURMAN
DMEK Preparation Forceps
for removal of Descemet layer with smooth round tips to prevent tearing the donor endothelium
overall length 75 mm



Geuder® 09 Forceps
DMEK, DSAEK, LASIK


▣-38628 YÖRÜK
DMEK Transplant Forceps
atraumatic grasping and removal of Descemet membrane
grasping area 6 mm
overall length 77 mm




▣-38634 GERTEN
DSAEK/DMEK Rhexis Forceps
for Descemetorhexis
reversed triangular tips to score and remove the DM from the host cornea
23 Gauge/ 0.6 mm curved shaft
overall length 142 mm



▣-33428 PFAEFFL
Flaporhexis Forceps
2.5 mm x 1.5 mm cross serrated tips, oval to lift Femto-LASIK flaps
overall length 103 mm



▣-31899 LASIK Corneal Flap Forceps
angled shafts, 9 mm
highly polished round and smooth shafts
overall length 108 mm




▣-18914 BURATTO
LASIK Corneal Flap Forceps
curved shafts
cross serrated tips
0.8 mm diameter
overall length 100 mm




Geuder® 09 Forceps
DMEK, DSAEK, LASIK

Geuder® 09 Forceps
Lenticule

♣-33960 RUPAL SHAH
Lenticule Forceps
3.5 mm serrated tips
20 Gauge/0.9 mm curved shaft
to retrieve lenticule during
SMILE procedures
6 mm titanium handle, blue-colored
overall length 144 mm




♣-33961 **Lenticule Forceps**
7.5 mm shafts
with fine grooves
angled to side
to retrieve lenticule during
SMILE procedures
overall length 99 mm




Geuder® 09 Forceps
Nucleus Prechoppers, Implantation




♣-31518 **Nucleus Spreading Forceps**
cross-action
3.0 mm textured paddle tips
for cracking the nucleus
overall length 112 mm




♣-31519 AKAHOSHI
Nucleus Spreading Forceps
Phaco Prechopper
4.5 mm delicate sharp tips, 1 mm profile
2.0 mm cross-action
for cracking, dividing, and
rotating the nucleus
overall length 105 mm




♣-31720 AKAHOSHI
Prechopper
for micro-incisions, straight shafts,
1.7 mm cross-action
round handle
overall length 125 mm



♣-33440 VON BARSEWISCH
Anterior Chamber Forceps
extra delicate
8.5 mm angled shafts
1.5 mm platform
overall length 72 mm




♣-32165 FECHNER
Implantation Forceps
smooth curved shafts
overall length 80 mm




Geuder® 09 Forceps
Nucleus Prechoppers, Implantation

GS-18935
DODEN
Implantation Forceps
curved
delicate lens holding shafts
overall length 110 mm



GS-31962
KOHLEN
Implantation Forceps
silicone tube coated
10 mm angled shafts
overall length 107 mm





Geuder® 09 Forceps
IOL

GS-18918
DICK
Fixation and Holding Forceps
for Artisan/Verisyse IOL
with guide mark for lens maximum hold
and stability
18 Gauge/ 12 mm curved shaft
overall length 143 mm



GS-18928
SEKUNDO
Implantation Forceps
for retropupillary implantation of
Iris-Claw-IOL
20 Gauge/ 0.9 mm curved shaft with
lens holding plate, three-point fixation
and IOL centration guide mark
6 mm titanium handle, blue-colored
overall length 143 mm




Geuder® 09 Forceps
IOL

GS-31902
Folder
foldable IOL insertion forceps
cross action
65 mm concave jaws
overall length 103 mm




GS-31903
SPALECK
Folder
foldable IOL insertion forceps
for silicone and acrylic IOLs
angled shafts, concave jaws
overall length 102 mm



GS-31908
Folder
IOL folding forceps
for easy pickup and folding acrylic IOLs
overall length 105 mm



GS-31994
H.-R. KOCH
Prefolder
for foldable silicone, hydrogel or
acrylic IOLs
overall length 135 mm




GS-31905
H.-R. KOCH
Folder
for foldable 3-piece IOLs
6.0 mm concave jaws
overall length 92 mm




Geuder® 09 Forceps
IOL


GS-31901 **Holding Forceps**
for foldable IOLs
highly polished round and smooth shafts
overall length 108 mm




GS-18941 **STROBEL Holding Forceps**
for Pharmacia CeeOn Edge IOL
overall length 110 mm



GS-32823 **Explantation Forceps**
for soft IOLs
2.0 mm x 1.4 mm paddle tip with
interlocking serrations and guide pin
overall length 108 mm




GS-32997 **SNYDER - OSHER IOL Forceps**
to hold the IOL while cutting it
inside the capsular bag
2 mm long jaws with
0.5 mm serrated tips
22 Gauge/ 0.7 mm lightly curved shaft
overall length 99 mm




Geuder® 09 Forceps
Superior Rectus


GS-18820 **REMKY Superior Rectus Forceps**
angled
1 x 2 teeth
overall length 109 mm




GS-18800 **ELSDCHNIG Superior Rectus Forceps**
standard, straight
1 x 2 teeth
overall length 110 mm



GS-19105 **LEYDHECKER Superior Rectus Forceps**
curved, tip width 0.8 mm
1 x 2 teeth, 1.1 mm
overall length 110 mm



GS-32335 **CASTROVIEJO Superior Rectus Forceps**
angled with tying platform
1 x 2 teeth, 0.9 mm
overall length 100 mm



Geuder® 09 Forceps

Muscle, Nasal, Cerclage, Retinal

GS-19112

**BANGERTER
Muscle Forceps**
10 teeth
self retaining
overall length 100 mm



GS-19120

GS-19125

**SATTLER
Muscle Forceps**
with slide lock
overall length 100 mm
3 teeth, bulbous tip, right
3 teeth, bulbous tip, left



GS-32350

**TROELTSCH
Nasal Forceps**
serrated tips
for packing and dressing during ENT procedures
overall length 150 mm




GS-19107

**ESSEN MODEL
Cerclage Forceps**
strongly curved
3 mm platform with dual grooves
overall length 104 mm



GS-33434

**PFENNIGWERTH
Holding Forceps**
for vitreoretinal lenses
curved with concave lens holding tips
overall length 110 mm




Geuder® 09 Forceps

Muscle, Nasal, Cerclage, Retinal

GS-33438

Forceps
for HEIMANN retinal tacks
overall length 100 mm




Geuder® 09 Forceps

Scleral Plug, Trocar


GS-33439

Scleral Plug Forceps
cross action handle, angled jaws
to securely hold plugs during insertion
and removal
overall length 105 mm



GS-33752

**WATZKE
Spreading Forceps**
to stretch silicone sleeves
cross action handle and
serrated tips
overall length 115 mm



GS-37430

Trocar Forceps
for trocar systems with groove
overall length 110 mm



Geuder® 09 Forceps
Vitreoretinal




☒-33831

REIMANN

Vitreoretinal Forceps

for foreign body removal
with retractable magnet and 3-prongs
14 Gauge/ 2.1 mm straight shaft
shaft length 35 mm
overall length 138 mm



Geuder® 09 Forceps
Vitreoretinal




Vitreoretinal Forceps

6 mm titanium handle
blue-colored
overall length 155 mm

straight
serrated tips


☒-36202

20 Gauge/ 0.9 mm



☒-36362
☒-36366


TANO
straight, serrated tips
23 Gauge/ 0.6 mm
25 Gauge/ 0.5 mm



straight
silicone coated



☒-36212

20 Gauge/ 0.9 mm



☒-36216
☒-36286


straight
grip tips
20 Gauge/ 0.9 mm
23 Gauge/ 0.6 mm



straight
serrated "crocodile"
tips



☒-36256
☒-36342
☒-36322

20 Gauge/ 0.9 mm
23 Gauge/ 0.6 mm
25 Gauge/ 0.5 mm



☒-36292
☒-36282
☒-36312


straight
endgripping,
extra delicate
0.2 mm x 0.05 mm
20 Gauge/ 0.9 mm
23 Gauge/ 0.6 mm
25 Gauge/ 0.5 mm



foreign body
forceps

☒-36262

20 Gauge/ 0.9 mm



Single-Use Interchangeable Forceps
available on page 176



Geuder® 09 Forceps
Chalazion, Ptosis

☒-18680
LAMBERT
Chalazion Forceps
round, 9 mm inside diameter
solid and fenestrated plates
overall length 92 mm



☒-18685
☒-18690
DESMARRES
Chalazion Forceps
oval
solid and fenestrated plates
20 mm, overall length 88 mm
25 mm, overall length 90 mm



☒-18700
CAUER
Chalazion Forceps
oval, 13 mm x 11 mm inside diameter
interlacing serrated tips
solid and fenestrated plates
overall length 99 mm



☒-18715
BERKE
Ptosis Forceps
22 mm long serrated tips
slide lock mechanism
overall length 110 mm




Geuder® 09 Forceps
Mosquito

☒-18170
☒-18171
HARTMANN
Mosquito Forceps
serrated
18 mm jaws
overall length 100 mm
straight
curved



☒-18180
☒-18181
HALSTEAD
Mosquito Forceps
serrated
25 mm jaws
overall length 125 mm
straight
curved



☒-18185
☒-18182
HALSTEAD
Mosquito Forceps
1 x 2 teeth
serrated
25 mm jaws
overall length 125 mm
straight
curved



Geuder® 09 Forceps
Towel, Serrefine



Ge-18160 BACKHAUS
Towel Forceps
overall length 90 mm



Ge-33200 SCHAEDEL
Towel Forceps
overall length 90 mm



Ge-18190 DIEFFENBACH
Serrefine
straight, overall length 35 mm



Ge-18195 DIEFFENBACH
Serrefine
curved, overall length 35 mm



Ge-31420 DIEFFENBACH
Serrefine
straight, overall length 50 mm



Ge-18197 DIEFFENBACH
Serrefine
extra delicate, overall length 50 mm

Geuder® 10 Forceps