

**INNOVATION**  
MADE BY GEUDER®

**NEW**

# PURE EFFICIENCY PHACO

Greatest Efficiency for  
Smallest Incisions

AUFFARTH  
MICS DOUBLE CROSS ACTION  
CAPSULORHEXIS FORCEPS



**G-32926 AUFFARTH  
MICS DOUBLE CROSS ACTION  
CAPSULORHEXIS FORCEPS**  
flat handle, cross-action,  
for incision size 1.6 mm,  
total length 100 mm



Prof. Gerd U. Auffarth, MD,  
Clinical Director  
University Eye Hospital Heidelberg

**More freedom of movement through the smallest clear cornea incision (down to 1.6 mm):**

- The narrow design guarantees great freedom of movement in the anterior chamber.
- The hand-crafted and precise forceps tips facilitate safe handling of the capsulorhexis.
- The overlying design of the forceps tips prevents the iris from getting caught in them, especially in cases of intraoperative floppy iris syndrome (IFIS).

**Geuder®**  
Precision made in Germany

GEUDER AG  
Hertzstrasse 4  
69126 Heidelberg  
Germany

Phone: +49 6221 3066  
Fax.: +49 6221 303122  
[info@geuder.de](mailto:info@geuder.de)  
[www.geuder.de](http://www.geuder.de)