

Medical Carriers User Guide



Tested to maintain temperature of 2-8 degrees

ICECATCH must be pre-conditioned
to -18 degrees before use

Product Name/Code	Carrier	ICECATCH	Placement	Ambient Temperature	*Temperature maintained for
Small Pathology Sample Carrier (PYTB1) Small Vaccine Carrier (PYTBVAC1)				 21°C 30°C	 17hrs 23min 12hrs 47min
Medium Pathology Sample Carrier (PYTB2) Medium Vaccine Carrier (PYTBVAC2)				 21°C 30°C 40°C	 29hrs 48min 14hrs 52min 8hrs 34min
Large Pathology Sample Carrier (PYTB3) Large Vaccine Carrier (PYTBVAC3)				 21°C 30°C	 31hrs 25min 15hrs 55min

*Study by TuV GmbH. ICECATCH preconditioned to -18°C for 24 hours. Ensure that placement of ICECATCH is correct with non-print insulated side in contact with contents. For more information please contact us.

