



### Long Handled Cheek Retractors

Excellent for photography.  
Autoclavable to 121°C. Sold in pairs.

Small	<b>DB04-0190</b>
Large	<b>DB04-0191</b>



### Clear Cheek Retractors



Very clear, excellent for photography.  
Cold sterilisable only. Sold in pairs.

Small	<b>DB04-0174</b>
Large	<b>DB04-0175</b>

### Opaque Cheek Retractors

Autoclavable up to 137°C. Sold in pairs.

Small	<b>DB04-0174PP</b>
Large	<b>DB04-0175PP</b>



### Kombee

Double ended occlusal and buccal cheek retractor. Keeps lips retracted for better occlusal pictures. Autoclavable up to 121°C. Sold as each.

Small	<b>DB04-0189</b>
Large	<b>DB04-0188</b>