

A COMPLETE

Neurotoxin Injection System



navaCLICK[™]
precision dosing

THE NAVACLICK™ INJECTION SYSTEM

Simplifies your preparation
& improves your efficiency



Connect the vial adapter for
needle-free reconstitution.



Attach the NavaClick for
needle-free withdrawal.



Secure the needle to
the NavaClick syringe.



Twist the plunger 90°
and click to inject.

Experience the *Brilliance*

Audible & tactile feedback for
EYES-ON-THE-PATIENT INJECTIONS

Metered-dose engineering for
PRECISION DOSING, EVERY CLICK

Low dead space needle & syringe for
REDUCED NEUROTOXIN WASTE

Sharp, ultra-thin needle for
ENHANCED PATIENT COMFORT

Sterile, self-sealing vial adapters for
NEEDLE-FREE DILUTION & DRAWS



“

“NavaClick is a game-changer for new injectors, and it makes an experienced injector even better.”

— Dr. Jason N. Pozner, MD

Precision dosing you won't have to see to believe

Allowing you to keep your eyes on the patient, the NavaClick announces – with an audible and tactile click – that you have delivered exactly .02 mL of neurotoxin.

“

“NavaClick is the best thing to happen to Botox since Botox.”

— Dr. Nayomi Omura, MD



Your patients will hear & feel the difference

Our trademark click and smooth ultra-sharp needles are gentle on the face and on the nerves, creating a new level of physical and emotional comfort for your patients. Don't be surprised when your patients — and their friends — request the syringe that clicks.



“

“I didn't think
injections could
be this easy.”


— Jenny, Patient

AVAILABLE NOW

The NavaClick™ Injection System



ORDER



LINEAGE | biomedical

lineagebiomedical.com