

the best in dental implant care




*Perfect choice
for your practice*



R138T
Montana Jack[®]
Titanium

PDT implant instruments feature lightweight, tactile-sensitive resin handles AND allow access without prosthesis removal for long-term success of implants.



strength

Our solid titanium won't scratch implants or leave debris behind -- and can be sharpened.



comfort

Reduce hand fatigue and pinch force while increasing tactile sensitivity with PDT's lightweight, wide handle design with knurling up to the shank.




efficiency

Our attention to blade size, angulation, curvature and balance allow for superior anatomical adaptation leading to overall comfort and efficiency.



WINGROVE™ GO-TO SET

QUICKLY ADAPT TO ANY IMPLANT




Your R900 Go-To Set includes:

- R661** **Wingrove™ B5-6 Ti**
- R680** **Wingrove™ L3-4 Ti**
- R693** **Wingrove™ N128B-L5 mini Ti**
- T216** **EasyView™ Probe UNC 12**
- T005A-F** **FlipTop™ Cassette**
- T061** **Diamond Head™
Ceramic Stone**

THE PINEYRO ARCH™ KIT

EXCELLENT HYGIENE ACCESS FOR
FULL FIXED IMPLANT RESTORATIONS



Your R930 KIT includes:

- R630** **PINEYRO ARCH™ TI 1**
Anterior Arch
- R631** **PINEYRO ARCH™ TI 2**
Anterior Arch
- R632** **PINEYRO ARCH™ TI 3**
Mesial/Distal Arch
- R633** **PINEYRO ARCH™ TI 4**
Specialty Arch




NO NEED FOR RESHARPENING



Montana Jack® Titanium R138T

TITANIUM AMAZING GRACEY™ MICRO MINIS

- Improved adaptation for small, narrow periodontal pockets
- 20% thinner and 1/2 the blade length of the Gracey Mini
- Extended terminal shank
- Strong enough to remove calculus and cement



GRACEY 1-2 MICRO MINI TI R694

GRACEY 11-12 MICRO MINI TI R695

GRACEY 13-14 MICRO MINI TI R696