





Top Facts about FORTUNA UNIVERSAL® OPTIFIX® HF:



- valve system without metal spring no contamination of dosing liquid
- suitable for analysis of micronutrients
- easy cleaning because of easy dismantling
- easy-to-use and solid volume setting system



- PTFE-encased piston with a sealing lip at its end for better tightness
- comes with one PTFE intake tube and two PTFE adaptors (GL 40 + 45)

Volume	grad.	Accuracy	Coefficent of Variation	ArtNo.
1 - 5 ml	0.2 ml	<u>≤ +</u> 0,6 %	<u>< +</u> 0.2 %	101.106-33
2 - 10 ml	0.2 ml	<u>≤ +</u> 0,6 %	<u>≤ +</u> 0.2 %	101.106-38
6 - 30 ml	1.0 ml	<u>≤ +</u> 0,6 %	<u>≤ +</u> 0.2 %	101.106-44
10 - 50 ml	1.0 ml	<u>≤ +</u> 0,6 %	<u>≤ +</u> 0.2 %	101.106-48



FORTUNA

UNIVERSAL