



Quantum VII BIA

SAFE, REPEATABLE, ACCURATE

The Quantum VII averaging stabilizes the effects of motion, respiration and circulation common in real time measurements.

Features:

Dual mode operation: Classical BIA or Scheduled BIA measurements

Applications include: assessing human body composition, motion, animals, plants, food and research

Averages and saves all BIA measurements with time stamp for spreadsheet analysis

Saved Classic and Scheduled BIA measurements are available from a built-in website using WIFI

Portable, handheld and light weight. Ideal for field BIA measurements

Long battery life on a single charge

Comes complete with everything needed for 100 tests

RJL Systems, Inc
33939 Harper Ave.
Clinton Township MI 48035 USA
800 528-4513
www.rjlsystems.com



Measuring and Managing Healthy Living