

# Lookin'Body<sup>120</sup> Spec Sheet

Get the most from your InBody. Use Lookin'Body.

Use Lookin'Body to manage all your InBody Test results. Lookin'Body connects to your InBody wirelessly and allows you to easily manage your InBody data on your personal computer. In addition, you can start InBody Tests remotely, track testing data, e-mail results to clients and customers, and share data on the same network.




## Convenient Features


**InBody Test** Lookin'Body allows you to control your InBody through your PC. Conduct the test and get results on your PC by clicking the InBody Test button.

**Data Management** With Lookin'Body, you can store and print results, e-mail result sheets, and edit client profiles. You can also customize reference ranges, result sheets, and logos for your clients. Lookin'Body allows you to quickly register a group of clients and export data as .csv or .jpeg for EMR. You can also add previous InBody Test data and view individual Health Reports.

**Health Report** The Health Report simplifies the InBody Test into clear, color-coded tables and graphs, making it easy to follow body composition trends line by line. It illustrates health history through accumulated measurements.




## Health Report



## Management

The screenshot shows the "Lookin'Body120" Management interface. It has tabs for "Select Member", "Select Test", and "Manage Results". The "Select Member" tab is active, showing a list of 44 members with columns for Name, ID, Height, Age, Gender, Member Info, and Health Report. A "User's Guide" section is visible on the right, stating "First select a member. \*Select Test or Manage Results on top after selecting." The interface also includes search and registration options.

## Results Sheet



## Member Info

The screenshot shows the "Member Info" interface. It includes a "Delete Member" button and a "Select Member" dropdown menu. The main area displays a form for "Required" information: Name (304800120), ID (304800120), Height (5 feet 10.5 inch), Gender (Male), Date of Birth (01/15/1977), and Age (37). There is also an "Optional" section for Mobile No., Telephone No., Address, Zip Code, E-mail, Registration Date (01.16.2014), and Memo. A "Save" button is at the bottom.