

SAM
MEDICAL

BLEEDING CONTROL KITS

BE PREPARED. PRESERVE LIFE.

Severe blood loss, often caused by motor vehicle accidents, acts of violence and natural disasters, can result in death within minutes. The SAM Medical Bleeding Control Kit includes essential medical devices designed to limit blood loss in emergency and traumatic situations. After calling emergency personnel, this life-preserving kit may help stabilize the injured in such events.



FEATURED TOURNIQUET

SAM XT EXTREMITY TOURNIQUET

Engineered for rapid application, SAM XT's innovative design normalizes the number of windlass turns needed compared to traditional tourniquets, enabling an easier, faster training of the product and intuitive use.






SAM[®] Bleeding Control Kit

KIT CONTENTS:

- 1 SAM XT EXTREMITY TOURNIQUET**
Slack is the main cause of failed tourniquet application. SAM XT's TRUFORCE™ Buckle technology auto-locks at a predetermined amount of circumferential force, eliminating nearly all tourniquet slack. Simply click, twist, and secure for a life-preserving tourniquet application.
- 2 MEDICAL GLOVES**
Nitrile protective medical gloves - Latex free, powder free, size L.
- 3 MEDICAL SHEARS**
Stainless steel medical shears allow the first responder to expose the wound and provide for appropriate treatment more efficiently.

- 4 COMPRESSED GAUZE**
Designed to treat wounds involving moderate to severe blood loss. Can be secured in place with Emergency Bandage. 100% cotton, 6-ply, sterile.
- 5 EMERGENCY BANDAGE**
Pressure bandage can be used directly over a wound or to hold gauze dressing in place with or without a pressure bar.
- 6 INSTRUCTION CARD**
Includes detailed instructions for use and a diagram for proper tourniquet placement.



PRODUCT SPECIFICATIONS

PACKAGED SIZE & WEIGHT*

Length: 9.7 in (24.5 cm)
 Height: 7 in (18 cm)
 Width: 3.2 in (8 cm)
 Weight: 10.9 oz (310 g)

OPERATING SPECIFICATIONS

After calling emergency personnel, this life-preserving kit may help stabilize the injured in traumatic events.

AVAILABLE MODELS

Vacuum packed - KT901-V-EN
 Non-Vacuum packed - KT901-NV-EN

GLOBAL SALES

U.S. Toll Free: 800.818.4726
 U.S. Local: 503.639.5474
 U.S. Fax: 503.639.5425
sammedical.com

	REF	KT901-V-EN
	REF	KT901-NV-EN
	SAM Medical Products [®] 27350 SW 95th Ave, Ste 3038 Wilsonville, OR 97070 USA	

* Size and weight vary by configuration and manufacturing process. ** Patent Pending
 SAM[®] is a registered trademark in the United States and other countries.
 ©2018 SAM Medical Products. All rights reserved.