

AED Plus Pelican Case

Total Protection for the First and Only Full-Rescue AED



Large Case
(19" x 15.4" x 7.6")
(48.3 cm x 39.1 cm x 19.3 cm)

Unbreakable, Watertight, Dustproof,
Chemical Resistant, and Corrosion Proof.

ZOLL
Advancing Resuscitation. Today.™

AED Plus Pelican Case

- Waterproof
- Crushproof
- Dustproof
- Corrosion proof
- Chemical resistant

Total Protection for the First and Only Full-Rescue AED

Providing an excellent shock-resistant carrier for your ZOLL® AED Plus®, the Pelican Case and its molded high-density foam provide total protection from the rigors of impact, vibration, and shock. Made of high-impact structural copolymer, it is extremely strong and durable. With its exclusive 0.25" (6.4 mm) neoprene O-ring, and ABS latches, this case seals perfectly. The Pelican Case also includes an automatic purge valve for quick equalization after changes in atmospheric pressure. The AED Plus Pelican Case is NATO-codified and tested to MIL C-4150J (Military Standard) and IP67 (water and dust ingress protection).

You break it, we replace it... forever.



A Rugged Case for the First and Only Full-Rescue AED

AED Plus Pelican Case Specifications

ZOLL offers two cases to meet your needs: a large case for easy storage of **CPR-D•padz™** and **pedi•padz® II** or a smaller case for the AED Plus alone. Each case comes complete with custom-molded foam inserts for a perfect fit.

Large Case (8000-0837-01)

Large Case with cutouts for AED Plus, **CPR-D•padz**, or **pedi•padz II**

Inside Dimensions: 17.9" x 12.8" x 6.8" (45.5 cm x 31.8 cm x 17.3 cm)

Outside Dimensions: 19" x 15.4" x 7.6" (48.3 cm x 39.1 cm x 19.3 cm)

Temperature Ratings: Min: -10° F (-23° C) Max: 210° F (99° C)

Buoyancy: Floats in salt water with 40 lbs. (18.1 kg) load

Color: Black

Approvals: IP67 Military Standard C-4150J

Small (8000-0836-01)

Small Case with cutouts for AED Plus only

Inside Dimensions: 14.7" x 10.4" x 6" (37.3 cm x 26.4 cm x 15.2 cm)

Outside Dimensions: 16" x 13" x 6.9" (40.6 cm x 33 cm x 17.5 cm)

Temperature Ratings: Min: -10° F (-23° C) Max: 210° F (99° C)

Buoyancy: Floats in salt water with 20 lbs. (9.07 kg) load

Color: Black

Approvals: IP67 Military Standard C-4150J

ZOLL Medical Corporation Worldwide Headquarters

269 Mill Road
Chelmsford, MA 01824
978-421-9655
800-348-9011

ZOLL Direct Sales Subsidiaries

Australia
+61 2 8424 8700
www.zoll.com.au

Austria
+43 650 413 6222
www.zollmedical.at

Canada
905 629 5005
www.zoll.com


France
+33 (1) 30 05 14 98
www.zoll.fr

Germany
+49 2236 87 87 0
www.zollmedical.de

The Netherlands
+31 (0) 488 41 11 83
www.zoll.nl

United Kingdom
+44 (0) 1925 846 400
www.zoll.com

For subsidiary addresses and fax numbers, as well as other global locations, please go to www.zoll.com/contacts.



ZOLL's AED Plus Large Pelican Case allows you to store your AED Plus and an extra set of **CPR-D•padz** with the option of also storing a set of **pedi•padz II**.

Specifications subject to change without notice.

©2007 ZOLL Medical Corporation. All rights reserved. "Advancing Resuscitation. Today." and **CPR-D•padz** are trademarks of ZOLL Medical Corporation. AED Plus, **pedi•padz II** and ZOLL are registered trademarks of ZOLL Medical Corporation.

Printed in U.S.A. 040710 9656-0173

ZOLL
Advancing Resuscitation. Today.™