


PRODUCT CATALOG

Prestan® Professional AED Trainer


Prestan® Professional AED Trainer


The **Prestan Professional AED Trainer** provides quality, durability and affordability like no other trainer in the market. Clear, confident voice instructions that provide easy-to-follow directions and preconfigured scenarios make training students easier. The Prestan Professional AED Trainer carries a three (3) year warranty.






- Includes 5 scenarios for varied training opportunities
- Automatic or Semi-automatic option
- Training Pads can be used 100 times during CPR/AED classes







- Module contains 2 languages
- Module contains the most up to date software and language capability in full compliance with the 2015 Guidelines









PRESTAN PROFESSIONAL AED TRAINER PARTS

PP-AEDT-101 (Single AED Trainer, English/Spanish) PP-AEDT-102 (Single AED Trainer, English/French) PP-AEDT-401 (4-Pack AED Trainers, English/Spanish)

PP-AEDT-402 (4-Pack AED Trainers, English/French) PP-APAD-1 (Single Set of Adult Pads)

PP-APAD-4 (4-Pack Sets of Adult Pads)

Optional Remote Control

- Volume, language and scenario control
- Simulate loose pad
- Pause, play or stop

• Initiate audible "Press Deeper" during CPR





Optional Pediatric Pads & Case






Prestan® AED ULTRATrainer™

Prestan® AED UltraTrainer™

The **Prestan AED UltraTrainer** provides the options instructors have been waiting for. Combining the most recent technology with realistic features not found on other AED Trainers, this product will excel in any training atmosphere. It is a compact and customizable unit providing quality, durability and affordability. With a three (3) year warranty, the **Prestan AED UltraTrainer** is sure to be an industry favorite.



- Pads are pre-connected and ready to use, unlike other AED trainers that use outdated plug-in technology
- Pads can be used 100 times (or more!) during CPR/AED classes
- One set of dual-graphics pads work for both adult and child training
- Pads have an embedded 'switch' which senses proper pad attachment to the manikin (Pad Sensing System) (patent pending)
- The silicone base adhesive on the pads, coupled with the foam structure, ensures easy application & removal with no residue



Prestan® Professional Training Manikins

Prestan® Professional Training Manikins

Realistic to the eye and the touch, the **Prestan Professional Training Manikin** series is unlike any other on the market. Available individually or in convenient multi-packs of four for class training purposes, our manikins are uniquely designed to give each student a realistic experience that provides for the best "real life" CPR outcomes. Prestan Manikins carry a three (3) year warranty.

MANIKIN SINGLES



Prestan Professional Adult Manikin PP-AM-100M-MS (monitored) PP-AM-100-MS (non-monitored)



Prestan Professional Child Manikin PP-CM-100M-MS (monitored) PP-CM-100-MS (non-monitored)



Prestan Professional Infant Manikin PP-IM-100M-MS (monitored) PP-IM-100-MS (non-monitored)

Manikin single units come with 10 Face-Shield/Lung-Bags and a nylon Carrying Case

- Available in Adult, Child and Infant models
- Realistic looking and performing manikin is more acceptable to both instructors and students.
- Far easier and faster to carry and setup at training sites
- Patented face/head tilt simulates how an actual victim's head would move if he required CPR and teaches students how to properly open the airway
- Easily inserted face-shield/ lung-bag provides visual chest rise to ensure students receive the most accurate and realistic training
- Face shield is secured so it won't slip and there's no need for adhesives
- Manikins are available with or without the CPR Feedback Monitor - Monitors may also be purchased separately to easily install into non-monitored Prestan Professional Manikins
- Prestan Manikins are not made with natural rubber latex
- Fully compliant with 2015 Guidelines

MANIKIN 4-PACKS



Prestan Professional Adult Manikin PP-AM-400M-MS (monitored) PP-AM-400-MS (non-monitored)



Prestan Professional Child Manikin PP-CM-400M-MS (monitored) PP-CM-400-MS (non-monitored)



Prestan Professional Infant Manikin PP-IM-400M-MS (monitored) PP-IM-400-MS (non-monitored)

Manikin 4-packs come with 50 Face-Shield/Lung-Bags and a nylon Carrying Case with a detachable Shoulder Strap



Prestan® Professional Training Manikins

AMAZING FEEDBACK TOOLS MAKE ALL THE DIFFERENCE!

Prestan's revolutionary *CPR Rate Monitor* provides instant feedback for the rate of chest compressions. Students can gauge speed on their own. Instructors can monitor several students quickly and easily. Lights continually signal rate progression, confirming compressions are within the 100-120 cpm range recommended by the AHA Guidelines. The yellow light will blink if the 120 cpm upper limit is exceeded, warning students to slow down.









Students can also be confident they are compressing to the appropriate depth. The unique *Clicker Mechanism* provides an audible confirmation when the correct compression depth is reached between 2.0 - 2.4 inches.

Clicker is in all Professional Manikins / CPR Rate Monitor is optional



Patented clamshell torso opening allows for quick and easy lung bag insertion and removal

ADULT JAW THRUST MANIKIN

Prestan's Jaw Thrust Head gives instructors the option to teach with the jaw thrust maneuver or head-tilt/chin-lift.





Jaw Thrust Heads may also be purchased separately - they are interchangeable with the heads of the Prestan Professional Adult Training Manikins

Prestan Professional **Adult Jaw Thrust**Training Manikin
PP-JTM-100M-MS (monitored)
PP-JTM-100-MS (non-monitored)
RPP-JTHEAD-1-MS (head only)

MANIKIN COLLECTIONS

Prestan Manikins are available in 2 convenient multi-packs.



The Prestan Collection (1 Adult / 1 Child / 1 Infant) PP-FM-300M-MS

Comes with 10 Face-Shield/Lung-Bags per Manikin, a nylon Carrying Case with a detachable Shoulder Strap and pre-installed Batteries for the Monitors



The Prestan Family Pack (2 Adults / 1 Child / 2 Infants) PP-FM-500M-MS

Comes with 10 Face-Shield/Lung-Bags per Manikin, a nylon Carrying Case with a detachable Shoulder Strap and pre-installed Batteries for the Monitors



Prestan Ultralite Manikin 4-Pack

Prestan® Ultralite® Manikins

Prestan is pleased to offer our newest and most portable manikin, the **Ultralite Manikin**. The Ultralite Manikin is amazingly easy to set up and is sold in a convenient and lightweight 4-pack for efficient training-on-the-go. This durable manikin offers an affordable means for CPR training, with all the quality and realism you expect from Prestan.











Prestan Ultralite Manikins PP-ULM-400-MS

- · Uniquely fast and easy set-up
- Stackable torsos and heads generate a compact package for transporting multiple units
- Clamshell opening provides easy access
- · Manikin is easy to clean
- Not made with natural rubber latex
- Compatible with AED Trainer pads
- · Realistic head tilt
- Visible chest rise with lung bag use
- Unique torso web design for realistic feel
- Anatomically correct design includes xyphoid process
- Realistic chest compressions to a 2" depth
- Life tested to a half-million compressions
- Three (3) year warranty

Ultralite Manikin 4-pack includes:

- 4 Manikin Torsos
- 4 Manikin Heads
- 4 Compression Pistons
- Nylon Piston Bag
- 50-pack Face-Shield Lung Bags
- · Nylon Carrying Case
- Detachable Shoulder Strap
- Instruction Sheet







Prestan® Products - Accessories

Prestan® Products - Accessories

DON'T RUN OUT!

Plan ahead to be sure you have enough Prestan accessories for your training needs. Prestan Accessories are available from your Prestan Authorized Distributor.

FOR YOUR MANIKINS



CPR Rate Monitor

for Professional Manikins CPR Feedback light system includes 120-upper-limit feature (2015 Guidelines)

> Adult: RPP-AMON-1 Child: RPP-CMON-1 Infant: RPP-IMON-1

Easily upgrade existing monitored manikins or install in manikins originally purchased without monitor



Face-Shield Lung Bags Adult 50-Pack: PP-ALB-50

Also available for Child, Infant, and **Ultralite Manikins**



Rescue Mask Adaptors

10-Bag: 10076-PPA 50-Bag: 10076-PPA-50



Face Shields

Adult 50-Pack: PP-AFS-50 Also available for Child, Infant, and **Ultralite Manikins**



Generic Face Shield

Adult 36-Count 10077-PPS

FOR YOUR AED TRAINERS



Replacement **Training Pads for the Professional AED Trainer** Single Set: PP-APAD-1



Replacement **Training Pads for the Professional AED Trainer**

4-Pack Set: PP-APAD-4



Replacement AED UltraTrainer

Training Pads for the Single Set: PP-UTPAD-1



Replacement **Training Pads for the AED UltraTrainer** 4-Pack Set: PP-UTPAD-4

Contact your Prestan Authorized Distributor for accessories and replacement parts (not all are listed above). Accessories are considered consumable items and are not covered by Prestan's three (3) year product warranty.





701 Beta Drive, Suite 3 • Mayfield Village, OH 44143 • 800-434-7989

Visit us at www.prestanproducts.com for more information and product resources

© Prestan Products LLC, January, 2018