

SAM Splint

SHAPE. SET. SECURE.

UNIVERSAL FLEXIBILITY

Touting unrivaled flexibility, SAM® Splint bends into any simple curve, becoming exponentially stronger and more supportive. SAM® Splint is globally acclaimed by emergency care providers, outdoor enthusiasts, the U.S. military, and even NASA.












ADAPTABLE C-CURVE™ PRINCIPLE

SAM® Splint leverages the revolutionary C-Curve™ Principle to mold its foam and aluminum layers into structural arcs, strong enough for pre- or post-cast care, while remaining lighter and more transportable than traditional splints that rely on heavy, rigid materials. Bend it back into its original form to stow away or remold for a separate injury.







SAM*Splint

ENGINEERED FOR VERSATILE SUPPORT

CUT-TO-SIZE

SAM® Splint's two layers of closed-cell foam and its aluminum interior can be cut with ordinary scissors to adapt to any size requirements.1

2 MOLDABLE ALUMINUM

Flat, O-temper aluminum inside the splint's outer layers molds easily, allowing for a wide spectrum of applications.

3 CLEANABLE EXTERIOR

Exterior closed-pore foam can be cleaned and disinfected with bleach and water for reuse.2

A RADIOLUCENT MATERIAL


Medical professionals are able to take X-rays and CT-scans without removing the splint from the fractured or injured limb.

WEATHER RESISTANT

Designed to function across the spectrum of global temperatures – from the Sahara to Mt. Everest.

6 LIGHTWEIGHT DESIGN

Substantially lighter than plaster or fiberglass splints, SAM® Splint requires only a wrap or tape to secure an injured bone or stabilize a joint.





PRODUCT SPECIFICATIONS

BLUE / ORANGE MODELS*

36" XL (SP1121FXL) 36" Roll (SP1121R)

36" Fold (SP1121F)

24" Industrial (SP1124F)

18" (SP1410F)

9" (SP1005F) Finger (SP1213L)

CHARCOAL MODELS*

36" XL (SP1109FXL) 36" Roll (SP1109R) 36" Fold (SP1109F)

18" (SP1412F)

9" (SP1107F)

OD GREEN MODELS*

36" Fold (SP1110F)

VET MODELS

36" VET (SP3003R)

GLOBAL SALES

U.S. Toll Free: 800.818.4726 U.S. Local: 503.639.5474 U.S. Fax: 503.639.5425


sammedical.com











SAM Medical Products® 27350 SW 95th Ave, Ste 3038 Wilsonville, OR 97070 USA

 $^{^{}f 1}$ Do not use a serrated blade. Ensure there are no sharp edges after cutting, and cover with tape.

² Care and maintenance guidelines are available at www.sammedical.com