

Which Prestan® AED Trainer is Right for Me?



> Professional vs. UltraTrainer™ <

In the market for a new AED Trainer? Which one best suits your training needs?



Prestan Professional AED Trainer

Features

Prestan AED UltraTrainer

Similar to Full-Size AED
Unit is 9.5" x 7.25" x 2.75"
and weighs 2.2 lbs.

Size and Weight

Compact and Portable
Unit is 7.25" x 5.25" x 2"
and weighs 1.2 lbs.

Yes

CPR Feedback

Yes

Yes

Pad Sensing System (PSS)

Yes

Yes

Pre-connected Pads

Yes

5 Scenarios Available

Shock Simulation

5 Scenarios Available

Yes

Auto / Semi-Automatic Shock

Yes

Yes

Rescue Breaths

Toggle On/Off

Yes

Compression Metronome

Toggle On/Off

No

Child Button

Yes

Available

Pediatric Pads

Dual Adult/Child Pads

Available

Remote Control

No

Included

Long-Life "C" Batteries

Included

Yes

Low Battery Prompt

Yes

3 Years

Prestan Warranty

3 Years

Visit www.prestanproducts.com for full details on these quality products.
Contact your Prestan Authorized Distributor for pricing and ordering.

