

PHILIPS

Emergency Care
and Resuscitation

HeartStart Intrepid



**Respond quickly.
Treat confidently.**

Philips HeartStart Intrepid

Options/accessories included with shipment

Part number	Description
867172	Base unit 3-lead ECG trunk cable (AAMI/IEC), 2.7 m 3-lead ECG leadset, snap AC power cord Lithium Ion battery Pads cable with test plug Therapy cable collar 1 roll paper 1 package electrodes Printed Quick Card User documentation CD Service manual CD

Monitoring options

These options cannot be combined

Option	Description
A01	SpO₂ monitoring Pulse oximetry Includes adult reusable fingertip transducer
A02	SpO₂, NBP monitoring Pulse oximetry, and non-invasive blood pressure Includes adult reusable fingertip transducer, adult blood pressure cuff, and NBP air hose
A03	SpO₂, NBP, EtCO₂ monitoring Pulse oximetry, non-invasive blood pressure, and capnography Includes adult reusable fingertip transducer, adult blood pressure cuff, NBP air hose, oral/nasal cannula, Capnotrak water filter, Capnotrak dehumidification tubing
A04	SpO₂, EtCO₂ monitoring Pulse oximetry, and capnography Includes adult reusable fingertip transducer, oral/nasal cannula, Capnotrak water filter, Capnotrak dehumidification tubing
A05	EtCO₂ monitoring Capnography Includes oral/nasal cannula, Capnotrak water filter, Capnotrak dehumidification tubing
A06	SpO₂, NBP, EtCO₂, Temp monitoring Pulse oximetry, non-invasive blood pressure, capnography and temperature Includes adult reusable fingertip transducer, adult blood pressure cuff, NBP air hose, oral/nasal cannula, Capnotrak water filter, Capnotrak dehumidification tubing

Note: Capnography consumables and temperature probes must be ordered separately

Therapy/diagnostic options

These options can be combined

Option	Description
B01	External pacing Adds noninvasive transcutaneous pacing
B02	AED mode Adds AED mode for adults and infant/child
B03	12-lead ECG acquisition Adds ability to acquire 12-lead ECG <i>Must be ordered with one of C13, C14, or C15</i>
B04	Data transmission Adds the ability to transmit data from the device
B08	Q-CPR CPR measurement and feedback Includes CPR meter 2, Q-CPR meter cable, Q-CPR meter pads

Accessory options

These options can be combined

Option	Description
C01	External paddles Adult and infant/child external paddles with PCI (patient contact indicator) Water resistant
C03	USB data drive
C04	Bedrail mount Includes one bedrail mounting bracket
C05	Lithium Ion battery Battery for use with HeartStart Intrepid only
C06	DC power module Medical DC - DC power supply
C20	Cable management straps Includes two cable straps
C21	Carry bag Intrepid accessory storage system Includes carry bag with side pouches, rear pouch and shoulder strap
C25	Docking station bundle Docking station, power cord, and bedrail hook

ECG accessory options

You may order only one of these options

Option	Description
C10	3-lead ECG, grabber 3-lead ECG trunk cable and lead set, grabber
C11	5-lead ECG, snap 5-lead ECG trunk cable and lead set, snap
C12	5-lead ECG, grabber 5-lead ECG trunk cable and lead set, grabber
C13	5/5 Snap ECG lead sets 10-lead trunk cable, 2 m 5-lead ECG, limb, snap 5-lead ECG, chest, snap
C14	5/5 Grabber ECG lead sets 10-lead trunk cable, 2 m 5-lead ECG, limb, grabber 5-lead ECG, chest, grabber
C15	3/7 Snap ECG lead sets 10 -lead trunk cable, 2 m 3-lead ECG, snap 7-lead ECG, snap

Data management options

Option	Description
D01	Intrepid wireless communication module Required for cellular/mobile data transmission

User documentation options

Option	Description
LP1	Instructions for use Includes a printed copy of the instructions for use

Service documentation options

You can only order a maximum of one of these options

Option	Description
SM1	Service manual, English Printed copy of the English service manual
SM2	Service manual, Simplified Chinese Printed copy of the Simplified Chinese service manual

Accessories and consumables

ECG monitoring

ECG monitoring electrodes

Part number	Description
M2202A	Adult radiotranslucent foam solid gel electrodes, 5 electrodes/pouch (300/case)
989803148821	Small Adult radiolucent foam, solid gel ECG electrode, 30 electrodes/pouch (600/case)

3-lead cables and lead sets

Part number	Description
M1669A	3-lead ECG trunk cable, AAMI/IEC
M1671A	3-lead ICU ECG set, grabber, AAMI
M1672A	3-lead ICU ECG set, grabber, IEC
M1673A	3-lead ICU ECG set, snap, AAMI
M1674A	3-lead ICU ECG set, snap, IEC
M1675A	3-lead OR ECG set, grabber, AAMI
M1678A	3-lead OR ECG set, grabber, IEC
989803173121	3-lead ECG set, disposable, bedside, AAMI
989803174201	3-lead ECG set, disposable, bedside, IEC

5-lead cables and lead sets

Part number	Description
M1602A	5-lead ICU ECG set, snap, chest, AAMI
M1604A	5-lead ICU ECG set, snap, chest, IEC
M1644A	5-lead ICU ECG set, snap, limb, AAMI
M1645A	5-lead ICU ECG set, snap, limb, IEC
M1668A	5-lead ECG trunk cable, AAMI/IEC
M1968A	5-lead ICU ECG set, grabber, limb, AAMI
M1971A	5-lead ICU ECG set, grabber, limb, IEC
M1973A	5-lead OR ECG set, grabber, limb, AAMI
M1974A	5-lead OR ECG set, grabber, limb, IEC
989803173131	5-lead ECG set, disposable, bedside, AAMI (1 pack=20 lead sets)
989803174211	5-lead ECG set, disposable, bedside, IEC (1 pack=20 lead sets)

ECG monitoring (continued)

12-lead cables and lead-sets

Part number	Description
M1663A	10-lead (5+5) ECG trunk cable, 2.0 m, AAMI/IEC
M1949A	10-lead (5+5) ECG trunk cable, 2.7 m, AAMI/IEC
M1976A	5-lead ICU ECG set, grabber, chest, AAMI
M1978A	5-lead ICU ECG set, grabber, chest, IEC
M1979A	5-lead OR ECG set, grabber, chest, AAMI
M1984A	5-lead OR ECG set, grabber, chest, IEC
M3526A	3-wire lead set plug, snap, AAMI
M3527A	7-lead ECG set, snap, chest, AAMI
M3528A	3-wire lead set plug, snap, IEC
M3529A	7-lead ECG set, snap, chest, IEC
989803176161	5-lead ruggedized ECG cable, snap, limb, AAMI
989803176171	5-lead ruggedized ECG cable, snap, chest, AAMI
989803176181	5-lead ruggedized ECG cable, snap, limb, IEC

Multifunction electrode pads

Part number	Description
M3501A	Adult/child multifunction defib pads, barrel connector, AAMI (10 sets)
M3504A	Infant multifunction defib pads, barrel connector, AAMI (5 sets)
M3713A	Adult/child plus pads (>10 kg) (10 sets)
M3716A	Adult/child radiolucent multifunction electrode pads (10 sets)
M3717A	Infant plus pads (<10 kg) (5 sets)
M3718A	Adult/child radiotransparent/reduced skin irritation multifunction pads (10 sets)
M3719A	Infant radiotransparent/reduced skin irritation multifunction pads (<10 kg) (5 sets)
989803139261	SMART Pads II for adults, children, and infants (1 set)
989803149981	SMART Pads III for adults, children, and infants (1 set)
989803149991	SMART Pads III for adults, children, and infants (5 sets)
989803158211	HeartStart adult multifunction electrode pads (1 set)
989803158221	HeartStart adult multifunction electrode pads (5 sets)
989803166021	HeartStart Adult/Child preconnect pads (10 sets)

External paddles

Part number	Description
M3543A	External paddles with PCI, water resistant
M4759A	Paddle electrode replacement for M3543A paddles
989803196431	Efficia external paddles with PCI, water resistant
989803197591	Paddle electrode replacement for 989803196431 paddles

Internal paddles

Part number	Description
M1741A	7.5 cm switchless internal paddles
M1742A	6.0 cm switchless internal paddles
M1743A	4.5 cm switchless internal paddles
M4741A	7.5 cm switched internal paddles
M4742A	6.0 cm switched internal paddles
M4743A	4.5 cm switched internal paddles
M4740A	Internal paddles adapter cable

Pads cables and adapters

Part number	Description
M3507A	Hands-free pads cable, barrel style
M3508A	Hands-free pads cable, plug style - 2.2 m (7 ft.)
989803197111	Pads adapter cable
989803158661	Hands-free pads cable, Q-CPR, HeartStart connector only, plug-style
05-10200	Pads adapter, converts barrel connector to plug-style

Paper

Part number	Description
989803138171	75 mm chemical thermal paper, 10 rolls

Data management

Part number	Description
989803202611	USB data drive
989803202921	Wireless communication module

Test loads and shorted plugs

Part number	Description
M3725A	Test load, 50 Ohm, plug-style
M1781A	Test load for use with M3507A pads cable, 50 Ohm, plug-style
989803171271	Test plug, shorted

Power

Part number	Description
989803202601	Lithium Ion battery
989803202931	Medical DC-DC power supply

SpO₂ sensors and cables

Part number	Description
M1131A	Single-patient SpO ₂ sensor, Pediatric/Adult finger (0.5 m)
M1132A	Single-patient SpO ₂ sensor, Infant finger (0.9 m)
M1133A	Single-patient SpO ₂ sensor, Neo/Infant thumb/toe, Adult finger (0.9 m)
M1134A	Single-patient SpO ₂ sensor, Neonatal foot/hand, Infant thumb/toe, Adult finger, adhesive-free (0.9 m)
M1140A	Single-patient Alar SpO ₂ sensor, Adult/Pediatric nasal
M1191B	Reusable SpO ₂ sensor, Adult finger (2.0 m)
M1191BL	Reusable SpO ₂ sensor, Adult (3.0 m)
M1192A	Reusable SpO ₂ sensor, Pediatric/Small Adult, glove (1.5 m)
M1194A	Reusable SpO ₂ sensor, Pediatric/Adult, ear clip (1.5 m)
M1196A	Reusable SpO ₂ sensor, Adult finger (3 m)
M1196S	Reusable SpO ₂ sensor, Adult clip (2 m)
M1941A	SpO ₂ extension cable, (2 m)
M1943A	SpO ₂ adapter cable (1 m)
M1943AL	SpO ₂ adapter cable (3 m)

Noninvasive blood pressure cuffs

NBP air hoses

Part number	Description
M1598B	NBP air hose, Adult/Pediatric (1.5 m)
M1599B	NBP air hose, Adult/Pediatric (3.0 m)

Multi-patient Comfort Care cuffs - single hose

Part number	Description
M1572A	Comfort Care cuff - Pediatric
M1573A	Comfort Care cuff - Small Adult
M1574A	Comfort Care cuff - Adult
M1574XL	Comfort Care cuff - Adult XL
M1575A	Comfort Care cuff - Large Adult

Noninvasive blood pressure cuffs (continued)

Multi-patient Easy Care cuffs - single hose

Part number	Description
M4552B	Easy Care cuff - Infant
M4553B	Easy Care cuff - Pediatric
M4554B	Easy Care cuff - Small Adult
M4555B	Easy Care cuff - Adult
M4556B	Easy Care cuff - Adult long
M4557B	Easy Care cuff - Large Adult
M4558B	Easy Care cuff - Large Adult X-long

Single care disposable cuffs

Part number	Description
989803182281	Single Care cuff - Pediatric (20 pack)
989803182291	Single Care cuff - Small Adult (20 pack)
989803182301	Single Care cuff - Adult (20 pack)
989803182311	Single Care cuff - Adult XL (20 pack)
989803182321	Single Care cuff - Adult large (20 pack)

EtCO₂

Part number	Description
989803198891	CO ₂ nasal cannula, large
989803198901	CO ₂ nasal cannula, medium
989803198911	CO ₂ nasal cannula, small
989803198921	CO ₂ /O ₂ nasal cannula, large
989803198931	CO ₂ /O ₂ nasal cannula, medium
989803198941	CO ₂ /O ₂ nasal cannula, small
989803198951	CO ₂ oral/nasal cannula, large
989803198961	CO ₂ oral/nasal cannula, medium
989803198971	CO ₂ /O ₂ oral/nasal, large
989803198981	CO ₂ /O ₂ oral/nasal, medium
989803198991	Airway adapter set, endotracheal diameter (ET) > 4.0 mm
989803199001	Airway adapter set, ET ≤ 4.0 mm
989803199011	Water filter assembly
989803199021	CapnoTrak CO ₂ sampling extension line, 60" (1.5 m)
989803199031	CapnoTrak O ₂ delivery extension line, 60" (1.5 m)
989803199041	CapnoTrak dehumidification tubing

Temperature

Part number	Description
21075A	Reusable esophageal/rectal temperature probe - Adult, pediatric
21076A	Reusable esophageal/rectal temperature probe - Infant
21078A	Reusable skin surface temperature probe - Adult, pediatric, infant
21082A	Disposable temperature probe extension cable (3.0 m)
21082B	Disposable temperature probe extension cable (1.5 m)
M1837A	Disposable esophageal/rectal temperature probe – 9 FR, 36" (20 pack)
21090A	Disposable esophageal/rectal temperature probe – 12 FR, 36" (20 pack)
21091A	Disposable sterile skin surface temperature probe, 31" (20 pack)
21093A	Disposable esophageal/stethoscope temperature probe – 12 FR, 20" (20 pack)
21094A	Disposable esophageal/stethoscope temperature probe – 18 FR, 20" (20 pack)
21095A	Disposable esophageal/stethoscope temperature probe – 24 FR, 20" (20 pack)
21096A	Disposable Foley catheter temperature probe – 16 FR, 17" (10 pack)

Temperature (continued)

Part number	Description
21097A	Disposable Foley catheter temperature probe – 18 FR, 17" (10 pack)
M2255A	Disposable Foley catheter temperature probe – 14 FR, 17" (10 pack)

Q-CPR accessories

Part number	Description
989803204811	Q-CPR meter 2
989803163291	Q-CPR adhesive pads (10 pack)

Cable management/accessory storage

Part number	Description
989803190351	Cable management straps
989803194281	Therapy cable collar
989803202941	Intrepid accessory storage system
989803202971	Intrepid rear accessory pouch
989803202981	Intrepid carry case shoulder strap

Mounting solutions

Part number	Description
989803190361	Bedrail/roll stand mount
989803203001	Intrepid docking station

Please note not all accessories and supplies are available in all countries.
The list of accessories is subject to change without notice.

Upgrades

- Each upgrade can only be ordered once per device
- For upgrades with options, you must order one of the available options

Part number	Upgrade	Options	Description
867288	SpO ₂ upgrade	D01 D02	Adds SpO ₂ if the device has no A options Adds SpO ₂ if the device has option A05 (EtCO ₂) <i>Other combinations are not possible</i>
867289	NBP upgrade	D01 D02	Adds NBP if the device has option A01 (SpO ₂) Adds NBP if the device has option A04 (SpO ₂ and EtCO ₂) <i>Other combinations are not possible</i>
867290	EtCO ₂ upgrade	D01 D02 D03	Adds EtCO ₂ if the device has no A options Adds EtCO ₂ if the device has option A01 (SpO ₂) Adds EtCO ₂ if the device has option A02 (SpO ₂ and NBP) <i>Other combinations are not possible</i>
867291	Temp upgrade	D01	Adds temperature if the device has A03 options (SpO ₂ , NBP, EtCO ₂) <i>Other combinations are not possible</i>
867292	AED mode upgrade	D01 D02	Adds AED mode if the device does not have B01 option (Pacing) Adds AED mode if the device has B01 option (Pacing)
867293	External pacing upgrade	D01 D02	Adds pacing mode if the device does not have B02 option (AED mode) Adds pacing mode if the device has B02 option (AED mode)
867294	12-lead acquisition upgrade		Adds 12-lead acquisition
867295	Q-CPR upgrade		Adds Q-CPR
867296	Cellular transmission upgrade		Adds cellular transmission
867297	Hands-free upgrade		Adds hands-free (pads) capability
867298	Data transmission upgrade		Adds data transmission capability

Note: The HeartStart Intrepid device serial number must be provided for all upgrades

