

UNIRAY | HD

by DigiMed



The **UNIVERSAL** Intraoral Sensor

One Size, For Pediatric and Adult Patients

Increase patient comfort and confidence with our 1.5 size sensor. Digital x-ray quality starts with high definition, and now with the UniRay HD, you can accommodate all of your patients, from pediatric to geriatric and everyone in between. The UniRay HD works with most major imaging software, including Dexis!

COMPATIBLE WITH



UniRay HD works in all leading imaging software and can be used in conjunction with existing sensors.



Your Practice Deserves To Purchase The Right Sensor!
 Let Our Team of Experts Help You Make the Best Choice

SENSOR TECHNICAL DATA

SENSOR TECH	CMOS
ACTIVE AREA	24 x 33mm
DEMENSIONS	29.2 X 39.5mm
BODY THICKNESS	4.8 mm
PIXEL SIZE	14.9 Microns
SCINTILLATOR	CSI



UNIRAY HD



- Proven durable engineering
- 1.5 size sensor
- Shock absorbing material
- Ergonomic design
- Integrates with Dexis & most branded imaging software

CORD



- 2M cable length
- Kevlar reinforced virtually unbreakable
- Replaceable
- Direct high speed 2.0 USB connectivity

UNIRAY HD KIT INCLUDES



Protective Sheaths

Rinn Positioning Kit

Wall Mount Holder

[DEALER TAG HERE]

Warranty
2 year warranty

Expedited Rental Service
Available for NO Downtime!

CONTACT YOUR DEALER TODAY!



Denterpriseintl.com