

QUICKRAY®
Trusted Since 1998




The
HASSLE FREE
Intraoral Sensor
Low Cost, Never Any Hidden Fees

The QuickRay Digital X-Ray Sensor series gives you excellent image clarity at the highest resolution with minimum radiation exposure; and connects directly to your computer's USB port without requiring a special interface box.



TWAIN compatible, QuickRay works in all leading imaging software and can be used in conjunction with existing sensors.

COMPATIBLE WITH






Your Practice Deserves To Purchase The Right Sensor!

Let Our Team of Experts Help You Make the Best Choice

SENSOR TECHNICAL DATA

SENSOR TECH	CMOS
RESOLUTION	Line Pairs >20 lp/mm
ACTIVE AREA	Size 1: 19.95mm x 30.02mm Size 2: 25.99mm x 36.00mm
SCINTILLATOR	CSI
PIXEL SIZE	19 x19 Microns
COMPATIBILITY	Windows 32/64 Bit & MAC OS
BODY THICKNESS	5.3 MM



QUICKRAY SENSOR



- Proven Durable Engineering
- Hermetically Sealed Housing
- Shock Absorbing Material
- Ergonomic Design
- Thin with Rounded Corners for Patient Comfort

CORD



- 2M Cable Length
- Kevlar Reinforced Virtually Unbreakable
- Replaceable
- Direct High Speed 2.0 USB Connectivity



ACTUAL IMAGE

SENSOR KIT INCLUDES



Protective Sheaths UNI-VERSE-ALL Positioning Kit Wall Mount Holder



Warranty

3 Year Warranty
7 Year Replacement Plan



Expedited Rental Service

Available for NO Downtime!

CONTACT US TODAY!



Outstanding Dental Equipment Since 1989

www.VideoDental.com | sales@videodental.com | 800.323.2690

156 East Granada Boulevard, Ormond Beach, FL 32176