

QUICKRAY^{DX}



The Affordable **HIGH RESOLUTION** Intraoral Sensor

Low Cost, Never Any Hidden Fees

The QuickRayDX Digital sensor offers superb image quality at a fraction of the cost of others. No need to sacrifice quality for price when you can have both with QuickRayDX.

COMPATIBLE WITH



QuickRay DX includes a built in TWAIN allowing for integration into all imaging software that is TWAIN compatible.





Your Practice Deserves To Purchase The Right Sensor!

Let Our Team of Experts Help You Make the Best Choice

SENSOR TECHNICAL DATA

SENSOR TECH	CMOS
RESOLUTION	Line Pairs >25 lp/mm
ACTIVE AREA	Size 1: 25.4 x 36.8 x 4.4 mm Size 2: 30.4 x 41.9 x 4.4 mm
SCINTILLATOR	CSI
PIXEL SIZE	1300x1800
COMPATIBILITY	Windows 32/64 Bit
DATA INTERFACE	USB 2.0 & 3.0



QUICKRAY^{DX} SENSOR



- Proven Durable Engineering
- Hermetically Sealed Housing
- Shock Absorbing Material
- Ergonomic Design
- Thin with Rounded Corners for Patient Comfort

CORD



- 2M Cable Length
- Kevlar Reinforced Virtually Unbreakable
- Replaceable
- Direct High Speed 2.0 USB Connectivity
- 70,000+ bendings RA test passed



ACTUAL IMAGE

SENSOR KIT INCLUDES



Protective Sheaths



Uni-Verse-All Positioning Kit



Wall Mount Holder



Warranty

2 Year Warranty
5 Year Replacement Plan



Expedited Rental Service

Available for NO Downtime!

CONTACT US TODAY!



VIDEODENTAL
CONCEPTS

Outstanding Dental Equipment Since 1989

www.VideoDental.com | sales@videodental.com | 800.323.2690
100 East Granada Boulevard Suite 219 Ormond Beach, FL 32176